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List of Experiments of Physics Laboratory – I

1. Determination of the Hall voltage developed across the sample material and also the Hall coefficient, mobility of charge carriers and carrier concentration of that material.
2. Determination of the resistivity and energy band gap of a semiconducting material using 4-Probe Method.
3. Determination of the excitation potential of Argon using Franck-Hertz apparatus.
4. Magnetic Susceptibility of Liquids – Quincke's Method, Powder-Faraday Method
5. Hysteresis (B – H Curve)
6. Solar-Cell Characteristics
7. Determination of Planck's Constant
8. Characteristics of Photo Diode, Photo Transistor, LDR, LED
9. To determine the wavelength of He- Ne LASER using Michelson interferometer
10. I-V characteristics of photoresistor.

MAULANA AZAD
NATIONAL INSTITUTE OF TECHNOLOGY, BHOPAL
DEPARTMENT OF PHYSICS

M.Sc. I Sem

EXPERIMENT: HALL EFFECT

AIM:

To determine the Hall voltage and Hall Coefficient of Germanium crystal.

APPARATUS

Hall Probe (Ge Crystal); Hall Effect Set-up, Electromagnet, Constant Current Power Supply, Digital Gaussmeter.

THEORY:

A static magnetic field has no effect on charges unless they are in motion. When the charges flow, a magnetic field directed perpendicular to the direction of flow produces a mutually perpendicular force on the charges. When this happens, electrons and holes will be separated by opposite forces. This means that the current of carriers will be deflected from a straight line in y -direction. In other words, there is a component of the velocity in y -direction and the surfaces perpendicular to the y -direction will become charged as soon as the current (or the magnetic field) is switched on. The flow-lines of the carriers will look like this,


Fig.1

They will in turn produce an electric field (E_h) which depends on the cross product of the magnetic intensity, H and the current density, J

$$E_h = RJ \times H \quad (1)$$

Where R is called the Hall Coefficient

Now, let us consider a bar of semiconductor, having dimension, x , y and z . Let J is directed along X and H along Z then E_h will be along Y . Then we could write

$$R = \frac{V_H / y}{JH}$$

Or

$$R = \frac{V_h}{I_x} \times \frac{z}{H} \quad (2)$$

Where V_h is the Hall voltage appearing between the two surfaces perpendicular to y and $I = I_{yx}$


Fig.2

The Hall Effect phenomenon is different for different charge carriers. In most common electrical applications, the conventional current is used partly because it makes no difference whether you consider positive or negative charge to be moving. But the Hall voltage has a different polarity for positive and negative charge carriers, and it has been used to study the details of conduction in semiconductors and other materials which show a combination of negative and positive charge carriers. n-type semiconductors are doped with trace quantities of a group V element (e.g. phosphorus, arsenic or antimony) as a pentavalent donor. p-type semiconductors are doped with trace quantities of a group III element (e.g. boron, aluminium or indium) as a trivalent acceptor. n-type semiconductors conduct by mobile negative charges (electrons), p-type semiconductors by mobile positive charges (holes). These electrons and holes are known as majority carriers. For semiconductors the Hall effect provides a means of determining the sign of the majority charge carriers and their number density.


Fig.3

FORMULA USED:

$$R = \frac{V_h}{I_x} \times \frac{z}{H}$$

Where V_h is Hall voltage, I_x current, H magnetic intensity and z is thickness of crystal.

$$R = \text{Slope} \times \frac{z}{H}$$

Slope ($\Delta V_h / \Delta I_x$) obtained from V_h versus I_x plot.

PROCEDURE:

- Connect Hall probe to set-up (Voltage and Current)
- Switch ON the setup and adjust current (say mA).
- Switch over the display to voltage side. See the Zero field potential.
- Place probe in the magnetic field as shown in fig. 3 and switch on electromagnet power supply and adjust the current of power supply at ~ 1.3 Amp (I_{em}) for 1000 Gauss constant magnetic field (measured by Gaussmeter).
- Vary current (I_x) of probe and measure the Hall voltage (V_h) corresponding to that current. Take 6 readings.
- Repeat experiment for different electromagnetic field (1100 Gauss and 1200 Gauss) by increasing the current (I_{em}) to electromagnet from power supply.
- Plot V_h versus I_x graph,
- Compute the slope ($\Delta V_h / \Delta I_x$).

OBSERVATIONS:

Thickness of probe (z) = 0.1 cm

OBSERVATION TABLE:

SET-1 Electromagnet power supply $H_1=1000$ Gauss ($I_{em} = \dots$ Amp)			SET-2 Electromagnet power supply $H_2=1100$ Gauss ($I_{em} = \dots$ Amp)			SET-3 Electromagnet power supply $H_3=1200$ Gauss ($I_{em} = \dots$ Amp)		
S.No.	I_x (mA)	V_h (mV)	S.No.	I_x (mA)	V_h (mV)	S.No.	I_x (mA)	V_h (mV)
1.			1.			1.		
2.			2.			2.		
3.			3.			3.		
4.			4.			4.		
5.			5.			5.		
6.			6.			6.		

CALCULATIONS:

Plot V_h versus I_x graph and compute slope ($\Delta V_h / \Delta I_x$). Calculate R for each set,

$$R_1 = \text{slope} \times \frac{z}{H_1}$$

$$= \dots\dots\dots \text{Volt.cm.amp}^{-1}.\text{G}^{-1}$$

Average,

$$R = \frac{R_1 + R_2 + R_3}{3}$$

$$= \dots\dots\dots \text{Volt.cm.amp}^{-1}.\text{G}^{-1}$$


Fig 4: Expected I-V curve

RESULT:

Hall Coefficient of Germanium crystal (R) = $\dots\dots\dots \text{Volt.cm.amp}^{-1}.\text{G}^{-1}$

PRECAUTIONS AND SOURCE OF ERROR:

1. The magnet power supply can furnish large currents at dangerous voltage levels; do not touch exposed magnet coil contacts.
2. Do not exceed magnet current of 4 A.
3. Do not exceed Hall probe current of 0.4 A

over 0794

Apparatus: SK012- FOUR PROBE
APPARATUS

Formula used: The energy band gap, E_g of a semiconductor is given by

$$E_g = \frac{2k \cdot 2.3026 \times \log_{10} \rho}{T^{-1}} \text{ eV}$$

Where, $k = 8.6 \times 10^{-5} \text{ eV/deg}$

ρ is resistivity of semiconductor sample

Given by $\rho = \rho_0 / f(W/s)$

Where, $\rho_0 = \frac{V}{I} \times 2\pi s$

W is thickness of sample,

s is the probe spacing,

$f(W/s)$ function is the correction factor,

V is the voltage across the two inner probes,

I is the current through the two outer probes.

Procedure:

1. The p-type Ge sample is factory fitted. If there is any need for replacing or mounting a sample then, proceed as per following procedure . Unscrew the pipe of the four probe arrangement and place the sample on the base plate of 4-probe arrangement and adjust the pipe such that the 4-point probe lie in the middle of the sample as shown in figure 3.


Fig. 3 shows four probe arrangement for mounting sample.

Gently applying some pressure and tighten the pipe in this position such that all the probes are in contact with sample.
(Warning! Applying excessive pressure may break the sample.)

2. Put the four probe arrangement in the oven and connect the lead of the oven to socket(10) as shown in figure 2. Also insert the PT-100 temp. sensor into the hole given at the top of the 4-probe arrangement.
 (Note: Do not switch 'ON' the oven. Keep the switch of oven(8) in OFF position.)
3. Connect the RED and BLACK plug leads of 4-probe arrangement to 4 mm sockets(7) marked as "Voltage".
4. Connect the YELLOW plug leads to 4mm sockets(6) marked as " Current".
5. Change switch (5) to current, so display (3) show current reading.
6. Switch 'ON' the apparatus. Slightly increase the current using current knob (2), say 4mA and note that the voltage should be positive. If it is not, interchange the current leads (6).
7. Set the current to desired value (say 8 mA) using current adjusting knob(2) as shown in figure 2. Also select the range of multiplier using switch (9) to X1 or X10 on voltage display. (It is always better to start with lower range i.e X1)
8. Switch 'ON' the oven using switch(8) as shown in figure 2. Green LED will glow ,showing that the oven is ON.
9. Change switch (5) to temperature, so display (3) show temperature.
10. Note the probe voltage on display (4) for different values of temperature as shown on display (3).

Observations:

- (1) Current, $I = 8.03 \text{ mA}$ (constant)
- (2) Distance between probes, $s = 0.24 \text{ cm}$
- (3) Thickness of sample, $W = 0.05 \text{ cm}$

S. No.	Temperature (°C)	Voltage (mV)	Temperature T(K)
1	23	119.1	296
2	30	124.3	303
3	40	131.2	313
4	50	138.2	323
5	60	145.2	333
6	70	151.4	343
7	80	155.5	353
8	90	155.0	363
9	100	147.5	373
10	110	130.8	383
11	120	106.1	393
12	130	81.3	403
13	140	62.9	413
14	150	47.9	423

S. No.	$T^{-1} \times 10^3$	ρ ($\Omega \text{ cm}$)	$\log_{10} \rho$
1	3.38	3.31	0.52
2	3.30	3.45	0.54
3	3.19	3.64	0.56
4	3.10	3.83	0.58
5	3.00	4.03	0.61
6	2.92	4.20	0.62
7	2.83	4.31	0.63
8	2.75	4.30	0.63
9	2.68	4.03	0.61
10.	2.61	3.63	0.56
11.	2.54	2.94	0.47
12.	2.48	2.26	0.35
13.	2.42	1.74	0.24
14.	2.36	1.33	0.12

3. Plot a graph ($T^{-1} \times 10^3$) versus $\log_{10} \rho$ and find the slope of the curve as shown in a typical graph figure 4


Fig.4 A typical graph showing the resistivity of Germanium crystal as a function of inverse temperature


So, the energy band gap of germanium is

$$E_g = 2k \cdot 2.3026 \times \log_{10} \rho$$

$$= 2k \times 2.3026 \times \frac{LM}{MN} \times 10$$

$$= 2 \times 8.6 \times 10^{-5} \times 2.3026 \times \frac{LM}{MN} \times 10$$

$$= 0.396 \times \frac{LM}{MN} \text{ eV}$$

From above graph the slope is found to be 1.78

Thus,

$$E_g = 0.396 \times 1.78 \text{ eV}$$

$$= 0.70 \text{ eV}$$

Result : Energy band gap for Ge is 0.70 eV

MAULANA AZAD
NATIONAL INSTITUTE OF TECHNOLOGY, BHOPAL
DEPARTMENT OF PHYSICS

EXPERIMENT: Franck-Hertz Experiment

Aim:

To determine the first excitation potential of gas (Argon) by Franck-Hertz experiment.

Apparatus:

Franck Hertz experiment- model FH3001, Tetrode tube filled with experimental Argon gas, filament.

Experimental Set-up:

The experimental set up involves a tube containing low pressure experimental gas fitted with four electrodes: an electron-emitting cathode (K), a mesh grid (G1) for minimizing space charge effects a mesh grid (G2) for acceleration, and an anode (A). The anode was held at a slightly negative electrical potential relative to the grid G2 (although positive compared to the cathode), so that electrons had to have at least a corresponding amount of kinetic energy to reach it after passing the grid and thereby making the dips in the plate current more prominent. Instruments were fitted to measure the current passing between the electrodes, and to adjust the potential difference (voltage) between the cathode (negative electrode) and the accelerating grids Fig (1).


Figure 1

Formula Used:

If V_n is the potential corresponding to n^{th} peak and V_1 is the potential corresponding to 1^{st} peak then

$$\text{Mean } 1^{\text{st}} \text{ excitation potential} = \frac{V_n - V_1}{n - 1}$$

Where $(n-1)$ is the number of dips between 1^{st} and n^{th} peak.

Theory

From the early spectroscopic work it is clear that atoms emit radiation at discrete frequencies; from Bohr's model, the frequency of the radiation ν is related to the change of energy levels through $E = h\nu$. It is then to be expected that transfer of energy to atomic electrons by any mechanism should always be in discrete amounts. One such mechanism of energy transfer is through inelastic scattering of low-energy electrons.

Franck and Hertz in 1914 set out to verify these considerations.

- (a) It is possible to excite atoms by low energy electron bombardment.
- (b) The energy transferred from electrons to the atoms always had discrete values.
- (c) The values so obtained for the energy levels were in agreement with spectroscopic results.

The Franck-Hertz experiment elegantly supports Niels Bohr's model of the atom, with electrons orbiting the nucleus with specific, discrete energies. Franck and Hertz were awarded the Nobel Prize in Physics in 1925 for this work.


Figure 2. Franck Hertz Experiment

Operating Principle:

The Frank-hertz tube in this instrument is a tetrode filled with the vapour of the experimental substance. Fig.1 indicates the basic scheme of experiment.

The electrons emitted by filament can be accelerated by the potential V_{G2K} between the cathode and the grid G_2 . The grid G_1 helps in minimizing space charge effects. The grids are wire mesh and allow the electrons to pass through. The plate (A) is maintained at a potential slightly negative with respect to the grid G_2 . This helps in making the dips in the plate current more prominent. In this experiment, the electron current is measured as a function of the voltage V_{G2K} . As Voltage increases, the electron energy goes up and so the electron can overcome the retarding potential V_{G2A} to reach the plate (A). This gives rise to a current in the ammeter, which initially increases. As the voltage further increases, the electron energy reaches the threshold value to excite the atom in its first allowed excited state. In doing so, the electrons lose energy and therefore the number of decreases. This decrease is proportional to the number of inelastic collisions that have occurred. When the V_{G2K} is increased further and reaches a value twice that of the first excitation potential, it is possible for an electron to excite an atom halfway between the grids, lose all its energy, and then gain a new enough energy to excite another atoms resulting in a second dip in the current. The advantage of this type of configuration of the potential is that the current dips are much more pronounced, and it is easy to obtain five-fold or even larger multiplicity in the excitation of the first level i.e. one can get 5 peaks (dips) or more.

The Franck-Hertz experiment elegantly supports Niel's Bohr's model of the atom, with electrons orbiting the nucleus with specific, discrete energies. Franck and Hertz were awarded the Nobel Prize in Physics in 1925 for this work.

Procedure:

1. Before the power is switched 'ON' make sure all the control knobs are at their minimum position and Current Multiplier knob at 10^{-7} or 10^{-8} or 10^{-9} (whichever suitable) position.
2. Switch 'ON' the power.
3. Turn the manual- Auto Switch to manual and check that the Scanning Voltage Knob is at its minimum position.
4. Turn Voltage Display Selector to V_{G1K} and adjust the V_{G1K} knob until voltmeter reads 1.5V.
5. Turn Voltage display selector to V_{G2A} and adjust the V_{G2A} knob until the voltmeter reads 7.5V.

When you have finished step 1-5, you are ready to do the experiment. Rotate V_{G2K} knob and observe the variation of plate current I_p with the increase of V_{G2K} . The current reading would show maxima and minima periodically. The magnitude of maxima could be adjusted suitably by adjusting the filament voltage and the value of Current Multiplier. Now take the systematic

(1/2 division). Plot the graph with output current I_p on Y-axis and accelerating voltage V_{G2K} at X-axis.

Observation Table :

V_{G1K} : 1.5V

V_{G2A} : 7.5V

S No.	Acceleration Potential V_{G2K} (Volts)	Plate Current I_p (nano Amperes)
1.		
2.		
3.		
..		
..		
..		



Results:

1. Graph with output current on Y-axis and accelerating voltage shows series of dips in current at approximately 12.1 volt (say) increments (fig 2).
2. At low potential differences—up to 12.1 volts when the tube contained argon vapour—the current through the tube increased steadily with increasing potential difference. The higher voltage increased the electric field in the tube and electrons were drawn more forcefully towards and through the accelerating grid.
3. At 12.1 volts the current drops sharply, almost back to zero.
4. The current increases steadily once again if the voltage is increased further, until 24.2 volts is reached (exactly 12.1+12.1 volts).
5. At 24.2 volts a similar sharp drop is observed.

Precautions:

1. During the experiment (manual), when the voltage is over 60V, please pay attention to the output current indicator, if the ammeter reading increase suddenly, decrease the voltage at once to avoid

- the damage of the tube.
- If you want to change the value of V_{G1K} , V_{G2A} and Filament Voltage during experiment, please first adjust the value of V_{G2K} to 'Zero'.
- Whenever the filament voltage is changed, please allow 2-3 minutes for its stabilisation.
- When the Frank-Hertz Tube is already in the socket, please make sure the following before the power is switched 'ON' or 'OFF', to avoid damage to the tube.
- Manual - Auto switch is on Manual and Scanning and Filament Voltage knob at its minimum position (rotate it anticlockwise) and current multiplier knob at 10^{-7} .

MAULANA AZAD
NATIONAL INSTITUTE OF TECHNOLOGY, BHOPAL
DEPARTMENT OF PHYSICS

M. Sc. I Sem

Quincke's Method

AIM:

To determine the magnetic susceptibility of Manganese sulphate solution by Quincke's method.

APPARATUS:

Quincke's tube with stand, Solution: MnSO_4 , Travelling Microscope, Electromagnet, Constant current Power Supply, Gaussmeter.

THEORY:

In electromagnetism, the magnetic susceptibility χ_m is a proportionality constant which is dimensionless and indicates the degree of magnetization of a material in response to an applied magnetic field. Mathematically, we express it as

$$\chi_m = \frac{M}{H} = \mu_r - 1 \dots \dots \dots (1)$$

Where M is the magnetization, H is the applied field and μ_r is the relative permeability of the material. The experimental set-up used for Quincke's method is shown below


Fig. Arrangement for Quincke's method

Here Manganese sulphate solution under investigation is placed in a vertical U-tube with one limb of wide bore and the other with narrow bore. The narrow limb is placed in between the pole pieces of the electromagnet. It should be noted that the surface of the liquid in the narrow limb must lie at the line of centres of the pole pieces when the field is off.

When the current is switched on a strong field is appeared at upper surface of the narrow column while the lower portion will be in a state of comparatively weak field. Hence a force will act upon the column and if the liquid is paramagnetic it will rise.

An expression for the magnetic susceptibility of the liquid as given below

$$\chi_{sol} = \frac{2gh(\rho - \sigma)}{\mu_0 H^2}$$

Where ρ = density of solution
 σ = density of air
 g = acceleration due to gravity
 h = rise in the surface of the liquid
 H = applied field
 μ_0 = permeability of vacuum

PROCEDURE:

- Put the tube on stand and fix it with clamp.
- Manganese sulphate solution is placed in a vertical U-tube.
- The narrow limb is placed in between the pole pieces of the electromagnet.
- Increase the value of current.
- Take the reading of magnetic field H from gaussmeter.
- Adjust the horizontal cross wire of the eye piece of microscope on the meniscus and note this reading of the microscope.
- Repeat the experiment with different value of current.
- Plot the graph between power supply current and magnetic field.
- Plot the graph between rise in liquid level h and H^2 .

OBSERVATIONS:

- Density of solution ρ = gm/cm³
- Density of air σ = gm/cm³
- Acceleration due to gravity g = 980 cm/sec²

$$\text{Least count of travelling microscope} = \frac{\text{Smallest division on main scale}}{\text{Total division on vernier scale}} = \dots\dots\dots \text{cm}$$

Table: for determining rising level of solution

S. No.	Power supply Current I (Amp)	Magnetic Field H (Tesla)	Magnetic Field H (K Gauss)	H^2 (K Gauss) ²	Initial position of meniscus h_1 (cm)	Final position of meniscus h_2 (cm)	Rise of solution $h = (h_2 - h_1)$ cm
1.							
2.							
3.							
4.							
5.							

CALCULATIONS:

$$\chi_{sol} = \frac{2gh(\rho - \sigma)}{\mu_0 H^2}$$

= _____

RESULT:

Magnetic susceptibility of the Manganese sulphate solution =

SOURCE OF ERROR AND PRECAUTIONS:

- i. Adjust the microscope to get clearer picture of meniscus.
- ii. Use lamp and magnifying glass while taking readings.
- iii. Let the liquid rise and become steady before taking readings of meniscus.
- iv. Reading should be taken carefully.

Aim: To study B-H curve and to find out the values of coercivity, retentivity and saturation magnetisation of experimental material. (commercial Nickel).

Apparatus Used: Set up for B-H curve, experimental material (commercial Nickel), CRO, connecting leads.


Figure 1

Theory:

Introduction: A precise knowledge of various magnetic parameters of ferromagnetic substances and the ability to determine them accurately are important aspects of magnetic studies. These not only have academic significance but are also indispensable for both the manufacturers and users of magnetic materials. The characteristics which are usually used to define the quality of the substance are coercivity, retentivity, saturation magnetisation and hysteresis loss. Furthermore, the understanding of the behaviour of these substances and improvement in their quality demand that the number of magnetic phases present in a system is also known. A B-H curve plots changes in a magnetic circuit's flux density as the magnetic field strength is gradually increased. The resulting shape indicates how the flux density increases due to the gradual alignment of the magnetic domains (atoms, that behave like tiny

to completely remove the remanence is called coercivity. If we continue to increase the negative magnetic field strength, the magnetic material will again reach saturation in the opposite direction, and the new curve will be a mirror image of the original curve. The complete B-H curve is then usually described as a hysteresis loop. The area contained within a hysteresis loop indicates the energy required to perform the 'magnetise-demagnetise' process. 'Soft' magnetic materials require relatively little energy to become magnetised and demagnetised and, so have 'narrow' hysteresis loops, whereas 'hard' magnetic materials require a great deal of energy and have 'wide' hysteresis loops.

So, B-H Curves and Hysteresis Loops are a valuable tool for comparing the characteristics and behaviour of different magnetic materials, in order to select them for an appropriate application.


Figure 3

DESIGN PRINCIPLE

When a cylindrical sample is placed coaxially in a periodically varying magnetic field (say by the solenoid) the magnetisation in the sample also undergoes a periodic variation. This variation can be picked up by a pickup coil which is placed coaxially with the sample. Normally, the pickup coil is wound near the central part of the sample so that the demagnetisation factors involved are ballistic rather than the magnetometric.

For the uniform field H_a produced, the effective field H acting in the cylindrical sample will be

$$H = H_a - NM \text{ where } M \text{ is the magnetisation, or}$$

$$H = H_a - \frac{NJ}{\mu_0}$$

(1)

where N is the normalised demagnetisation factor including 4π and J is the magnetic polarization defined by

$$B = \mu_0 H + J \quad (2)$$

with $B = \mu H$ or $\mu_0(H + M)$ as magnetic induction. The signal corresponding to the applied field, H_a , can be written as

$$e_1 = C_1 H_a \quad (3)$$

where C_1 is a constant.

Further the flux linked with the pickup coil of area A_c due to sample of area A_s will be

$$\phi = \mu_0(A_c - A_s)H' + A_s B$$

Here H' is the magnetic field, in the free from sample area of the pickup coil, will be different from H and the difference will be determined by the magnitude of demagnetising field. However, when the ratio of length of the sample rod to the diameter of the pickup coil is more than 10, the difference between H and H' is too small, so that

$$\begin{aligned} \phi &= \mu_0(A_c - A_s)H + A_s B \\ &= \mu_0 A_c H + A_s(B - \mu_0 H) \end{aligned}$$

$$\Rightarrow \phi = \mu_0 A_c H + A_s J \quad (4)$$

The signal e_2 induced in the pickup coil will be proportional to $\frac{d\phi}{dt}$

After integration the signal (e_3) will, therefore be

$$e_3 = C_3 \phi = C_3 \mu_0 A_c H + C_3 A_s J \quad (5)$$

Solving equations (1), (3) and (5) for J and H give

$$C_1 C_3 A_c \left(\frac{A_s}{A_c} - N \right) J = C_1 e_1 - \mu_0 C_3 A_c e_1 \quad (6)$$

$$\text{and } C_1 C_3 A_c \left(\frac{A_s}{A_c} - N \right) H = C_3 A_s e_1 - \frac{N C_1 e_1}{\mu_0} \quad (7)$$

Based on these equations an electronic circuit may be designed to give the values of J and H and hence the hysteresis loop.

In case the sample contains a number of magnetically different constituents, the loop obtained will be the algebraic sum of individual loops of different phases. The separation of these is not

easy in a J - H loop while in a second derivative of J , $\frac{d^2J}{dt^2}$ the identification can be made very clear.

EXPERIMENTAL DESIGN AND ANALYSIS

The aim is to produce electrical signals corresponding to J and H as defined in Eqs. (6) and (7) so that they can be displayed on CRO (cathode ray oscilloscope).

Moreover, it should be able to display $\frac{d\phi}{dt}$ & $\frac{d^2\phi}{dt^2}$ as a function of H or usual time

base of the CRO.

The magnetic field has been obtained with a multilayered solenoid driven by the AC mains at 60 Hz and supplied through a variable transformer arrangement. The magnetic field has been calibrated with a Hall probe and is found to be within $\pm 3\%$ of the maximum value over a length of 5 cm. across the central region. The instantaneous current producing the field is passed through a resistor R_1 in series with the solenoid and measured with an AC ammeter. The resulting signal e_1 is applied across a 500Ω helipot and an adder amplifier through a $100\text{ K}\Omega$ resistance.

The signal e_2 corresponding to the rate of change of flux is obtained from a pickup coil wound on a non-conducting tube. Necessary arrangements have been made

COMMONLY USED PHRASES

to place the sample coaxially with the pickup winding and in uniform magnetic field. The pick-up coil is connected to point B (Fig. 2) through twisted wires, where e_2 constitutes the input for further circuit. To obtain J , e_2 is fed to an adjustable gain integrator. Because of capacitive coupling of pickup coil and solenoid, self inductance of pickup coil and integration operation an additional phase will be introduced in the output signal e_3 , whose sign can be made negative with respect to e_1 by interchanging the ends of the pickup coil. To render e_3 completely out of phase with e_1 , a phase shifter consisting of a $1\text{K}\Omega$ potentiometer and $1\mu\text{F}$ capacitor has been connected at the output of integrator. Amplitude attenuation due to this network is compensated by the gain of the integrator and is not important as addition of signals is performed afterwards.

The out of phase signals e_1 and e_3 are added at the input of a unity gain adder amplifier and its output which is proportional to J is applied to Y-input of a CRO. Fractions of these signals corresponding to the demagnetisation factor and area ratio form the input of another adder amplifier with gain 10 whose output after further amplification of 10 is fed to the X-input of CRO and gives H . It may be mentioned that the gains of the amplifier can be adjusted but should always be such that the operational amplifiers are not loaded to saturation.

The selector switch (SW) can change the Y-input of CRO to J , $\frac{dJ}{dt}$ or $\frac{d^2J}{dt^2}$. The $\frac{dJ}{dt}$ signal is taken directly from the pickup while $\frac{d^2J}{dt^2}$ is obtained through an operational amplifier differentiator.

Let us now analyze the circuit.

The magnetic field at the centre of the solenoid for current i flowing through it will be

$$H_0 = Ki \quad (8)$$

$$\text{also } e_1 = R_1 i \quad (9)$$

with symbols defined above Eq. (9) reduces to Eq. (3) with $C_1 = R_1/K$. Further, when the sample is placed in a pickup coil of n turns

$$\begin{aligned} e_2 &= n \left(\frac{d\phi}{dt} \right) \\ &= n\mu_0 A_c \left(\frac{dH}{dt} \right) + nA_s \left(\frac{dJ}{dt} \right) \end{aligned} \quad (10)$$

by substituting ϕ from Eq. (4), we get

$$\begin{aligned} -e_3 &= -g_1 \int e_2 dt \\ &= -g_1 n \mu_0 A_c H - g_1 n A_s J \end{aligned} \quad (11)$$

Where g_1 is the gain of the integrator and phase shifter combination. The sum of e_1 and $-e_3$ after amplification becomes.

$$\begin{aligned} e_y &= -g_y (e_1 - e_3) \\ &= -g_y (C_1 H - g_1 n \mu_0 A_c H + C_1 \frac{NJ}{\mu_0} - g_1 n A_s J) \end{aligned} \quad (12)$$

Using Eq. (1), (3) and (11), g_y is the gain of this amplifier. If we adjust $C_1 = g_1 n \mu_0 A_c$, then

$$e_y = g_y g_1 n A_c \left(\frac{A_s}{A_c} - N \right) J \quad (13)$$

Fraction α and β of e_1 and $-e_3$ respectively, are added together at the input of the first amplifier for the X-input. If g_x be the total gain of both amplifiers we get

$$\begin{aligned} e_x &= g_x (e_1 - \beta e_3) \\ &= g_x g_1 n \mu_0 A_c (-\beta) H + g_x g_1 n A_c (N - \beta \frac{A_s}{A_c}) J \end{aligned} \quad (14)$$

after substituting $C_1 = g_1 n \mu_0 A_c$, J will be eliminated from the right hand side of (14). By adjusting α and β such that

$$\alpha = \frac{A_s}{A_c} \text{ and } \beta = N \quad (15)$$

$$\text{we get } e_x = g_x g_1 n \mu_0 A_c \left(\frac{A_s}{A_c} - N \right) H \quad (16)$$

Equation (13) and (16) can be written as

$$H = G_0 \frac{e_x}{\left(\frac{A_s}{A_c} - N \right)} \quad (17)$$

$$\text{and } J = \frac{\mu_0 g_x e_y}{g_y \left(\frac{A_s}{A_c} - N \right)} \quad (18)$$

Where

$$\frac{1}{G_0} = g_x g_1 n \mu_0 \Lambda_c \quad (19)$$

Equations (17) and (18) define the magnetic quantities H and J in terms of electrical signals e_x and e_y respectively.

METHOD

Calibration

When an empty pickup coil is placed in the solenoid field, the signal e_2 will only be due to the flux linking with coil area. In this case $J = 0$, $\alpha = 1$, $N = 0$ so that $H = H_a$ and Eqs. (13) and (16) yield

$$e_y = 0 \quad \text{and} \quad e_x = G_0^{-1} H_a \quad (20)$$

i.e. on CRO it will be only a horizontal straight line representing the magnetic field H_a . This situation will, obviously, be obtained only when the condition for (13) is satisfied. Thus without a sample in the pickup coil a good horizontal straight line is a proof of complete cancellation of signals at the input of the Y-amplifier. This can be achieved by adjusting the gain of the integrator and also the phase with the help of network meant for this purpose. From known values of H_a and the corresponding magnitude of e_x , we can determine G_0 and hence calibrate the instrument. The dimensions of a given sample define the values of demagnetisation factor and the area ratio pertaining to the pickup coil. The demagnetisation factor can be obtained from the Appendix. These values are adjusted with the value of 10 turn helipot provided for this purpose. The value of the area ratio can be adjusted upto three decimal places whereas that of N upto four (Zero to 0.1 max.). The sample is now placed in the pickup coil. The plots of J , $\frac{dJ}{dt}$ and $\frac{d^2J}{dt^2}$ against H can be studied by putting the selector switch at appropriate positions. The graph of these quantities can also be obtained from time base by using the internal time base of CRO.

Since eddy currents are present in conducting ferromagnetic materials, the resulting J-H loop has a small loop in the saturation portion due to difference in phases for the forward paths. Moreover, these plots do not show horizontal lines at saturation and hence their shapes can't be employed as a criterion to adjust the values of demagnetisation factor.

The values of loop width, intercept on the J-axis and saturation position are determined in terms of volts for different applied fields. Plots of these against magnetic field are then used to extract the value of coercivity, retentivity and saturation magnetic polarization. The first corresponds to the intercept of the width against currents straight line on the Y-axis and it is essentially the measure of the width under no shielding effects. On the other hand, the remaining two parameters are derived from asymptotic extensions of the corresponding plots because these refer to the situation when shielding effects are insignificant. Caution is necessary in making the straight line fit for loop widths as a function of current data as the points for small values of magnetic current have some what lower magnitudes. This is due to the fact that incomplete saturation produces lower coercivity values in the material. The geometrically obtained values of potentials are, in turn, used to find the corresponding magnetic parameters through equations (17) and (18).

If the area ratio for a particular sample is so small that the loop does not exhibit observable width, the signal e_x can be enhanced by multiplying α and β by a suitable factor and adjusting the two helipots accordingly. The ultimate value of the intercept can be normalised by the same factor to give the correct value of coercivity.

Observations

For this equipment diameter of pickup coil = 3.21mm

$$g_x = 100$$

$$g_y = 1$$

Sample : Commercial Nickel

Length of sample : 39 mm

Diameter of sample : 1.17 mm

Therefore,

$$\text{Area ratio } \left(\frac{A_s}{A_c} \right) = 0.133$$

Demagnetisation factor (N) = 0.0029 (Appendix)

Calibration

Settings : Without sample. Oscilloscope at D.C. Time base EXT. H Bal., Phase and DC Bal. adjusted for horizontal straight line in the centre. Demagnetisation at zero and Area ratio 0.40 at magnetic field 200gauss (rms)

$$e_x = 64\text{mm, or}$$

$$= 7.0\text{V (if read by applying on Y input of CRO)}$$

For Area ratio 1

$$e_x = 160\text{mm, or}$$

$$= 17.5\text{V}$$

From Eq. (20)

$$G_0(\text{rms}) = \frac{200}{160} = 1.25\text{gauss/mm}$$

$$G_0(\text{peak to peak}) = 1.25 \times 2.82 \\ = 3.53\text{gauss/mm,}$$

also

$$G_0(\text{rms}) = \frac{200}{17.5} = 11.43\text{ gauss/volt}$$

$$G_0(\text{peak to peak}) = 11.43 \times 2.82 \\ = 32.23\text{gauss/volt}$$

By adjusting N and $\frac{A_s}{A_c}$ as given above the J-H loop width is too small. Thus both are adjusted to three times i.e. 0.399 and 0.0087 respectively (full value of area ratio pot. = 1.000 and full value of demag. pot. = 0.100)

(a) Coercivity

S.No.	Mag. Field (rms) (Gauss)	2x Loop width (mm)
1.	30	7.0
2.	62	9.0
3.	94	11.0
4.	138	12.5
5.	179	14.0
6.	226	15.5
7.	266	16.75
8.	302	18.0
9.	336	18.75

(b) Saturation magnetisation

S.No.	Mag. Field (rms) (Gauss)	Tip to tip height (mv)
1.	29	205
2.	61	370
3.	96	400
4.	137	420
5.	176	430
6.	223	440
7.	264	445
8.	298	450
9.	331	450

(c) Retentivity

S.No.	Mag. Field (rms) (Gauss)	2xIntercept (mV)
1.	29	170
2.	61	260
3.	95	265
4.	136	270
5.	175	270
6.	219	275
7.	263	275
8.	302	275
9.	335	275

From the graphs Fig. (4) and (5)

Loop width = 2.9mm (after dividing by the multiplying factor 3)

2xIntercept = 280mV

Tip to tip height = 457.5mV

CALCULATIONS

(a) Coercivity

Since $e_x = \frac{1}{2} \times \text{loop width} = \frac{1}{2} \times 2.9 = 1.45 \text{ mm}$

$$H = \frac{G_0 e_x}{\left(\frac{A_s}{A_c} - N\right)} = \frac{3.53 \times 1.45}{(0.133 - 0.0029)} = 39.30e \text{ from equation (17)}$$

(b) Saturation magnetisation

$$\mu_s = \frac{J_s}{4\pi}$$

due to equation (2)

$(e_y)_s = \frac{1}{2} \times \text{tip to tip height} = 457.5/2 = 228.75 \text{ mV}$

$$\mu_s = \frac{J_s}{4\pi} = \frac{G_0 \mu_0 g_x (e_y)_s}{g_y \left(\frac{A_s}{A_c} - N\right) \times 4\pi}$$

from equation (18)

$$= \frac{32.23 \times 1 \times 100 \times 0.229}{1 \times (0.133 - 0.0029) \times 12.56} = 452 \text{ gauss}$$

(c) Retentivity

$$\mu_r = \frac{J_r}{4\pi}$$

due to equation (2)

$(e_y)_r = \frac{1}{2} \times (2 \times \text{Intercept}) = \frac{1}{2} \times 280 = 140 \text{ mV}$

$$\mu_r = \frac{J_r}{4\pi} = \frac{G_0 \mu_0 g_x (e_y)_r}{g_y \left(\frac{A_s}{A_c} - N\right) \times 4\pi} = \frac{32.23 \times 1 \times 100 \times 0.140}{1 \times (0.133 - 0.0029) \times 12.56} = 276 \text{ gauss}$$

Note : The above observation and calculation are given as a typical example. Test results of individual unit are supplied with the unit separately

MAULANA AZAD
NATIONAL INSTITUTE OF TECHNOLOGY, BHOPAL
DEPARTMENT OF PHYSICS

M.Sc. I Sem

EXPERIMENT: SOLAR CELL

AIM:

To plot the V-I Characteristics of the solar cell.

APPARATUS:

Solar cell, voltmeter, ammeter, light bulb, connecting wires etc.

THEORY:

The solar cell is a semi conductor device, which converts the solar energy into electrical energy. It is also called a photovoltaic cell. A solar panel consists of numbers of solar cells connected in series or parallel. The number of solar cell connected in a series generates the desired output voltage and connected in parallel generates the desired output current. The conversion of sunlight (Solar Energy) into electric energy takes place only when the light is falling on the cells of the solar panel. Therefore in most practical applications, the solar panels are used to charge the lead acid or Nickel-Cadmium batteries. In the sunlight, the solar panel charges the battery and also supplies the power to the load directly. When there is no sunlight, the charged battery supplies the required power to the load.


Fig. 1a Working principle of a solar

Fig. 1

A Solar cell operates in somewhat the same manner as other junction photo detectors. A built-in depletion region is generated in that without an applied reverse bias and photon adequate energy create hole-electrons pairs. In the solar cell, as shown in Fig. 1, the pair must diffuse a considerable distance to reach the narrow depletion region to be drawn out as useful current. Hence, there is higher probability of recombination. The current generated by separated pairs increases the depletion region voltage (Photovoltaic effect). When a load is connected across the cell, the potential causes the photocurrent to flow through the load.

The e.m.f. generated by the photo-voltaic cell in the open circuit, i.e. when no current is drawn from it is denoted by V_{oc} (V-open circuit). This is the maximum value of

e.m.f.. When a high resistance is introduced in the external circuit a small current flows through it and the voltage decreases. The voltage goes on falling and the current goes on increasing as the resistance in the external circuit is reduced. When the resistance is reduced to zero the current rises to its maximum value known as saturation current and is denoted as I_{sc} , the voltage becomes zero. A V-I characteristic of a photovoltaic cell is shown in Fig. 3.

The product of open circuit voltage V_{oc} and short circuit current I_{sc} is known as ideal power.

$\therefore$ Ideal Power = $V_{oc} \times I_{sc}$
The maximum useful power is the area of the largest rectangle that can be formed under the V-I curve.


Fig.2


Fig.3

PROCEDURE:

- Keep the solar cell in sunlight for 15 or 20 minutes so that it gets activated.
- Complete the circuit connections as shown in the figure (2) below. Then place the light 'S' at a distance of 15 cm from solar cell.
- Now find out the open circuit voltage V_{oc} by opening the connecting wire joining -ve of milli-ammeter to the load (i.e. connecting +ve of cell to the +ve of voltmeter and -ve of the cell to the -ve of voltmeter).
- Once again complete the connections as shown in fig. 2.
- Bring load selector switch at S.C. position and note the reading in the ampere meter under short circuit condition (if meter shows out of scale then decrease the light intensity). This maximum current is called short circuit current I_{sc} .
- Now introduce the load resistance in the circuit (start from low value of resistance) and note down the current and voltage reading.

- g) Repeat the same with different values of the load resistance.
 h) Draw a graph between V and I (fig. 3).
 i) Repeat the experiment with other intensities placing the light source at different positions.

OBSERVATIONS:

- a) Voltmeter reading for open circuit, $V_{oc} = \dots$ Volts
 b) Milliammeter reading with zero resistance, $I_{sc} = \dots$ mA.

Table: for determining voltage and current

S.No.	Voltage	Current density (J) mA/cm ²	Power density (mW/cm ²)
1.			
2.			
3.			
4.			
5.			

CALCULATIONS:

Value of $V_m = \dots$ volts

Value of $I_m = \dots$ mA

Maximum useful power = $V_m \times I_m$ mW

Ideal power $V_{oc} \times I_{sc} = \dots$ mW

RESULT:

The I-V characteristics of the solar cell is plotted on the graph paper using experimental readings.

SOURCE OF ERROR AND PRECAUTIONS:

- The solar cell should be exposed to sun light before using it in the experiment.
- Light from the lamp should fall normally on the cell.
- A resistance in the cell circuit should be introduced so that the current does not exceed the safe operating limit.
- Make sure that connections are made properly and ensure good contact.

MAULANA AZAD
NATIONAL INSTITUTE OF TECHNOLOGY, BHOPAL
DEPARTMENT OF PHYSICS

M. Sc. I Sem.

EXPERIMENT: PLANCK'S CONSTANT

OBJECT:

To determine the Planck's Constant using LED.

APPARATUS USED:

Variable Voltage source (0-2V DC), Current meter (0-20 mA/2000μA),
Temperature controlled Oven (ambient to 60°C), LED (Red, Yellow).

THEORY:

The basic idea in this measurement is that the photon energy, which from Einstein's relation is $E_f = h\nu$ is equal to the energy gap E_g between the valance and conduction bands of the diode. Energy gap is in turn equal to the height of energy barrier eV_o . Those electrons have to overcome to go from the n-doped side of diode junction to the p-doped side when no external voltage V is applied to the diode. In the p-diode side they recombine with holes releasing the energy E_g as photon with $E_f = E_g = eV_o$. Thus a measurement of V_o indirectly yields E_f and planck's constant. If ν is known or measured. However there are practical and conceptual problems in actual measurement.

Let us consider the LED diode equation:

$$I \propto \exp(-V_o / V_r) [\exp(V / V_r) - 1], \quad V = V_m - RI \quad (1)$$

$$\text{where } V_r = \frac{\eta k T}{e}$$

K = boltzmann constant, T = absolute temperature & e = electronic charge.
 V_m is voltmeter reading in external diode circuit and R is the contact resistance.
The constant η is material constant, which depends on type of diode, location of recombination region etc. The energy barrier eV_o is equal to the gap energy E_g when no external voltage V is applied. The quantities, which are constant in an LED, are impurity atom density, the charge diffusion properties and the effective diode area. The 'one' in the rectifier is negligible if $I \geq 2nA$, and the equation becomes,

$$\begin{aligned} I &\propto \exp[(V - V_o) / V_r] \\ &\propto \exp[e(V - V_o) / \eta k T] \end{aligned} \quad (2)$$

The height of potential barrier is obtained by directly measuring the dependency of diode current on the temperature keeping the applied voltage and thus the height of barrier fixed. The external voltage is kept fixed at a value lower than the barrier.

In our experimental set-up the variation of current I with temperature is measured over about a range of about 30°C at a fixed voltage $V(=1.8\text{volts})$ kept slightly below V_o . The slope of $\ln I$ vs $1/T$ curve gives

$$V_o = V - \left[\frac{\Delta \ln I}{\Delta T^{-1}} \times \frac{K}{e} \times \eta \right]$$

The constant η may be determined separately from I-V characteristic of the diode at room temperature from the relation

$$\eta = (e/kT)(\Delta V / \Delta \ln I) \quad (3)$$

(NOTE: value of η is already given in observation table)

The Planck's constant is then obtained by relation

$$h = eV_o \lambda / c \quad (4)$$

The contact resistance of LED is usually around 10ohm , while overall internal resistance of LED at applied voltage (1.8V) is few hundred ohms. The factor RI in expression $V=V_m-RI$ may therefore be neglected.

FORMULA:

The Planck's Constant,

$$h = \frac{e \times V_o \times \lambda}{c}$$

Where,

$$V_o = V - \left[\frac{\Delta \ln I}{\Delta T^{-1}} \times \frac{K}{e} \times \eta \right]$$

slope of curve of graph, $\frac{\Delta \ln I}{\Delta T^{-1}}$

K = Boltzmann constant

e = Electronic charge

η = Material constant

λ = Wavelength of LED

c = velocity of light

PROCEDURE:

Dependence of current (I) on temperature (T) at constant applied voltage to find Temperature coefficient

1. Keep the mode switch on V-I side and adjust the voltage across LED slightly below the band-gap of LED ie 1.8V for both yellow and red.
2. Change the mode of two-way switch to T-I side.
3. Insert LED in the oven and connect the other end of LED in the socket provided on set up. Before connecting the oven check that oven switch is in OFF position and SET Temperature knob is at minimum position. Now 1st display would read ambient temperature.
4. Set the different temperature 20, 25, 30,..... °C with the help of set –Temperature knob. Allow about 5 minutes on each set temperature to stabilize and take the reading of temperature and current.
5. Find the inverse of temperature and draw the graph between in I & (1/T).

OBSERVATIONS:

- Sample: (RED / Yellow) LED
- Voltage = 1.803V (constant for whole set of readings)
- Material constant of Red LED (η) = 1.755
- Material constant of Yellow LED (η) = 1.782
- Boltzmann Constant $K = 1.38 \times 10^{-23}$ J/Kelvin
- Electric Charge $e = 1.602 \times 10^{-19}$ Coulombs
- Wavelength of yellow light $\lambda = 5800 \times 10^{-8}$ cm
- Wavelength of Red light $\lambda = 6400 \times 10^{-8}$ cm

Table: Determination of temperature coefficient of Current

S.No.	Temperature (°C)	Temperature (°K)	1/T (K ⁻¹)	Current I (mA)	lnI (I in mA)
1.	30				
2.	35				
3.	40				
4.	45				
5.	50				
6.	55				
7.	60				

Draw graph between $\ln I \rightarrow 1/T$

CALCULATIONS:

Slope of curve


$$\frac{\Delta \ln I}{\Delta T^{-1}} =$$

$$V_o = V - \left[\frac{\Delta \ln I}{\Delta T^{-1}} \times \frac{K}{e} \times \eta \right]$$

The Planck's Constant,

$$h = \frac{e \times V_o \times \lambda}{c}$$

$$h = \dots \dots \dots \text{Joule.sec}$$

RESULT:

The Planck's Constant $h = \dots \dots \dots$ Joule.sec

PRECAUTIONS:

- i. V-I characteristic of LED should be drawn at very low current upto = 1000 μ A only, so that disturbance to V_o is minimum.
- ii. In T-I mode, make sure that the oven switch is 'OFF' and SET temp knob is at minimum position before connecting the oven.
- iii. On each setting of temperature, please allow sufficient time for the temperature to stabilized, between 5-6 minutes.

EXPERIMENT: Photodiode, Phototransistor, LDR, LED

AIM:

To study IV Characteristics of Photo Diode, phototransistor, LDR, and LED.

APPARATUS:

Photo diode, Photo transistor, LDR, and LED Kit

THEORY:

Photovoltaic and Photoconductive devices like Photodiodes, Photo resistors (also known as Light Dependent Resistors or LDRs), Phototransistors etc. are example of semiconductor Light Sensors. A Light Sensors or Photo Sensors is a device that measures the intensity of light, usually the Electromagnetic Radiation in wavelength range from Ultraviolet to far infrared.


PHOTODIODE

It is a form of light sensor that converts light energy into electrical energy (voltage or current). Photodiode is a type of semi conducting device with PN junction. Between the p (positive) and n (negative) layers, an intrinsic layer is present. The photo diode accepts light energy as input to generate electric current.

Internally, a photodiode has optical filters, built in lens and a surface area. When surface area of photodiode increases, it results in less response time. Few photo diodes will look like Light Emitting Diode (LED). It has two terminals as shown below. The smaller terminal acts as cathode and longer terminal acts as anode.

Figure 1

PHOTODIODE



Figure 2

The symbol of the photodiode is similar to that of an LED but the arrows point inwards as opposed to outwards in the LED as shown in Figure 2.

PHOTOTRANSISTOR

A Phototransistor is a Photo junction device that is similar to a regular transistor except that it has a light sensitive Base terminal (or Collector – Base Junction). In other words, a photo transistor can be considered as a Photodiode with current amplifier. A phototransistor



converts photons to charge directly, just like a photodiode, and in addition to this, a phototransistor also provides a current gain. The symbol of phototransistor is shown below.

The principle of operation of a phototransistor is similar to a photodiode in combination with an amplifying transistor. The light incident on the base of a phototransistor will induce a small current. This current is then amplified by normal transistor action, which results in a significantly large. Usually, when comparing to a similar photodiode, a phototransistor can provide a current that is 50 to 100 times that of a photodiode.

Characteristics of Phototransistor

As phototransistors are basically Bipolar NPN Transistors with large base – collector junction, the characteristics of a Phototransistor are similar to that of a simple BJT. Phototransistors are available as two-leaded or three-leaded devices. In a two-leaded Phototransistor, the Base terminal electrically unavailable and the device is completely dependent on light. Collector terminal is usually at a higher potential than Emitter to induce reverse bias at base – collector junction. When there is no light falling on the phototransistor, a small amount of leakage current known as Dark Current flows from collector to emitter. When there is enough light falling on the base terminal, a base current is produced, which is proportional to the intensity of the light.

LDR

A Light Detector or a Light Sensor is a device or circuit that detects the intensity of the light incident on it. When the light falls on the LDR, the light stays off and when the light stops falling on LDR, the LED glows.



When subjected to light energy, a Photoconductive light sensor will change its physical property. Photo Resistor is a common type of photoconductive device. Photo resistor is a semiconductor device that uses light energy to control the flow of electrons and therefore the flow of current in them.

The most common type of photoconductive cell is a Light Dependent Resistor or LDR. As the name implies, a Light Dependent Resistor is a semiconductor device that

changes its electrical resistance depending on the presence of light. A Light Dependent Resistor changes its electrical resistance from a high value of several thousand Ohms in the dark to only a few hundreds of Ohms when light is incident on it by creating electron - hole pairs in the material.

The most common material used to make a Light Dependent Resistor is Cadmium Sulphide (CdS). Other materials like Lead Sulphide (PbS), Indium Antimonide (InSb) or Lead Selenide (PbSe) can also be used as the semiconductor substrate.


Cadmium Sulphide is deposited as a thread pattern on an insulator in the shape of a zigzag line as shown below.


The reason for zigzag path is to increase dark resistance and therefore decrease the dark current. This cell is encapsulated in a glass to protect the substrate from contamination.

LED

Light Emitting Diode or simply LED is one of the most commonly used sources of light nowadays. Unlike (almost) legacy filament bulbs, LEDs (and fluorescent bulbs) need a special circuit to make them work. They are simply called as LED Drivers. The two most significant semiconductor light emitting sources extensively used in various applications are LASER diodes and LED's. The principle operation of LASER diodes is based on stimulated emission, whereas LED is based on spontaneous emission. A light emitting diode consists of two principal elements of semiconductor. They are positively charged P-type holes and negatively charged N-type electrons.

When the positive P side of the diode is connected to power supply and N side to the ground, then the connection is said to be in forward bias, which allows the electric current to flow through the diode. The majority and minority charge carriers in the depletion layer at the PN junction. The migration of electrons and holes in turn releases some amount of photons, which discharges energy in the form of monochromatic light at a constant wavelength. The color spectrum of LED emission is typically extremely narrow. The selection of emission of color from the LED is fairly limited due to the nature of semiconductor used in the manufacture. Commonly available colors of LED are red, green, blue, yellow, amber and white. The working voltage of red, green, amber and yellow colors is around 1.8 volts. The actual range of working voltage of a light emitting diode can be determined by the breakdown voltage of semiconductor material involve in the construction of LED. The color of the light emitted in LED is determined by the semiconductor materials that form the diode's PN junction.

PROCEDURE:

Phototransistor:

1. Focus the Light on the photo transistor, take reading of Volts and micro ammeter.
2. Change the Volt and distance of the light source and again take other set of reading and draw the graph.

RESULT:

The IV Characteristics of Photodiode, Phototransistor, LDR, and LED are shown in the graph.

SOURCE OF ERROR AND PRECAUTIONS:

1. Make sure that all the corrections are right and then switch ON the power.
2. Do not provide light continuously to the Photo Diode, phototransistor, LDR, and LED for more than 30 Mins.
3. Handle the apparatus carefully.

MAULANA AZAD
NATIONAL INSTITUTE OF TECHNOLOGY, BHOPAL
DEPARTMENT OF PHYSICS

B.Tech. I Sem/II sem

EXPERIMENT: MICHELSON INTERFEROMETERS

OBJECT:

To determine the wavelength of He-Ne laser with the help of Michelson interferometer

APPARATUS USED:

Michelson interferometer, laser source, screen etc.

THEORY:

The Michelson interferometer operates on the principle of division of amplitude rather than on division of wave front.

The important parts of a Michelson interferometer include a sturdy base, a diffusing glass, a beam splitter, a movable mirror with a micrometer screw for measuring distance of movement, a fixed mirror, and compensating glass. These parts are shown in a typical physical arrangement in Figure 1.


Fig.1 A typical Michelson interferometer

The light source shown to the left can be a white light or spectral light source. It is not generally considered a part of the Michelson interferometer.

An accurately machined micrometer screw is attached to the movable mirror, permitting the mirror to be moved toward or away from the beam splitter by a precisely determined amount. The micrometer generally has about one inch of movement and usually can measure mirror motion to within 0.002 mm. The two mirrors, beam splitters, and compensating glass all are made of carefully annealed glass (e.g., borosilicate-crown).

The optical arrangement for the Michelson interferometer is shown schematically in Figure 2. Light from a *broad* spectrum source *S*, such as an incandescent bulb, passes through a diffusing glass *DG* (e.g., a ground glass plate), and strikes the beam splitter *P*. The beam splitter is half-silvered glass plates (silvered on the back side) which reflects half of the light toward mirror *M*₁ and transmits half of the light (but the entire cross section) toward mirror *M*₂.

If the two mirrors are precisely aligned such that planes are exactly perpendicular to one another, thus ensuring that differences over different regions of the mirrors are constant, the fringe pattern will be seen by the observer at consist of a series of concentric rings Figure 3.

The fringes shown are called fringes of equal inclination. When the mirror *M*₁ is moved so as approach the condition for zero path difference, the fringe pattern will appear collapse with all fringes moving toward the center, then disappearing. Collapse of one fringe at the center will be equivalent to a displacement of mirror *M*₁ by an amount of λ . Thus if *N* is the number of fringes collapsing at the center when *M*₁ is moved by a distance of *d*₀ then

$$2d_0 = N\lambda$$

Using this relation λ can be obtained.

FORMULA USED:

Wavelength λ ,

$$\lambda = 2 d_0 / N$$

where *d*₀ is distance measured with micrometer screw and *N* is number of fringes collapsing

PROCEDURE:

- Set up the Michelson interferometer
- Mount the laser on the optical table with the beam parallel to the table surface
- The moveable mirror will be the one in the direct path of the laser beam.
- Align the mirrors using the laser beam.


Fig.2 Optical arrangement and light path in Michelson interferometer


Fig. 3 Circular fringes (equal inclination) seen in Michelson interferometer

- Adjust the distance of M_1 and M_2 from P to be nearly equal.
- Adjust the center of circular fringes in the field of view.
- Note down the least count of micrometer screw provided with the mirror M_1 .
- Record the initial reading of the screw. Now, move it solely and count the number of fringes collapsing at the center in steps of 25 and each time note down the micrometer reading in the tabular form.
- Measure the wavelength of the He-Ne laser using eq.3 .

OBSERVATIONS:

Least count of interferometer Main scale = cm

Least count of circular scale V_1 = cm

Least count of vernier scale V_2 = cm

TABLE:

No. of fringes collapsed	Main scale (cm) mm	Circular scale V_1 (cm) mm	Vernier scale V_2 (cm) mm	Total (cm) (mm)	Separation d_o for $N=50$ (cm) 25 mm
0					
25					
50					
75					
100					

Mean value of d_o = (cm)

CALCULATIONS:

Wavelength λ ,

$$\lambda = 2 d_o / N$$

$$= \text{..... cm}$$

$$= \text{..... } \text{\AA}$$

RESULT:

The wavelength of given He-Ne laser is $\text{\AA}$

PRECAUTIONS AND SOURCE OF ERROR:

- The two mirrors should be highly silvered on their front surfaces.
- The glass plate should be of same thickness.
- Adjustment should be made by moving M_1 only



EXPERIMENT I-V CHARACTERISTICS

Aim : To determine the current – voltage(I-V) characteristics of given Conductor and Semiconductor

Apparatus Required : Apparatus required A diode, A DC voltage supplier, Bread board, 100 ohm resistor, 2 multimeter for measuring current and voltage, connecting Wires

Theory

The increase (by Ohm's precise relationship between voltage and current) is different for different components and can be shown on an I-V graph :

If an electric conductor is connected to a voltage source with voltage V a current I is produced. We define resistance being the ratio of the voltage applied and current flowing through conductor:

$$R = V/I$$

The resistance depends on the material and geometry of the conductor. In the simplest case the current is proportional to the voltage applied. Thus, in this case the resistance is constant and the conductor termed Ohmic resistor (or conductor). For such resistors, Ohm's law holds true:

$$R = V/I \quad \text{at constant Temperature } T.$$

The I-V graph for a metallic conductor at constant temperature e.g. resistor, is very simple:

- The current is directly proportional to the potential difference.
- This is demonstrated by the straight line graph through the origin.

The I-V graph for a semiconductor diode is slightly different. A diode is used in a circuit to allow current to flow only in a specific direction:

- When the current is in the direction of the arrowhead symbol, this is forward bias. This is shown by the sharp increase in potential difference and current on the right side of the graph

- When the diode is switched around, it does not conduct and is called reverse bias. This is shown by a zero reading of current or potential difference on the left side of the graph.

OBSERVATION TABLE

For Conductor:

S.NO.	Voltage (V)	Current (mA)
1.		
2.		
3.		
4.		
5.		

For Semiconductor:

S.NO.	Voltage (V)	Current (mA)
1.		
2.		
3.		
4.		
5.		

Result : For I-V characteristics, the graph will be drawn between current and voltage.

List of Experiments of Physics Laboratory – II

1. To calculate the frequency of oscillations of RC phase shift oscillator
2. Study of RC Coupled CE amplifier – Two stages with feedback – Frequency response and voltage gain
3. Study of Push-pull amplifier using complementary – symmetry transistors power gain and frequency response.
4. Study of Active filters – low pass and high pass-first and second order frequency response and roll off rate.
5. Study of gain of Inverting/ Non-inverting amplifier and also study the frequency response characteristics and find out the bandwidth..
6. Study of frequency response characteristics of differentiator/ integrator.
7. Study of frequency response characteristics of Clipper Clamper circuit.
8. Verification of truth table of OR, AND and NOT gates.
9. To determine the frequency and wavelength in a rectangular waveguide working in TE₁₀ mode
10. To study the drain, transfer; drain resistance, amplification factor, and Trans conductance characteristics of an FET.

EXPERIMENT: PHASE SHIFT OSCILLATOR

AIM:

To determine the frequency of Phase Shift Oscillator.

APPARATUS:

Transistor, Capacitor, Resistor, CRO, Connecting wires.


THEORY:

An oscillator is an electronic circuit for generating an AC signal voltage with a DC supply as the only input requirement. The frequency of the generated signal is decided by the circuit elements used. An oscillator requires an amplifier, a frequency selective network and a positive feedback from the output to the input.

RC-Phase shift Oscillator has a CE amplifier followed by three sections of RC phase shift feedback Networks. The output of the last stage is returned to the input of the amplifier. The values of R and C are chosen such that the phase shift of each RC section is 60° . Thus The RC ladder network produces a total phase shift of 180° between its input and output voltage for the given frequencies. Since CE Amplifier produces 180° phase shift, the total phase shift from the base of the transistor around the circuit and back to the base will be exactly 360° or 0° . This satisfies the Barkhausen condition for sustaining oscillations. Total loop gain of this circuit is greater than or equal to 1, this condition used to generate the sinusoidal oscillations.

The frequency of oscillations of RC-Phase Shift Oscillator is,

$$\frac{1}{2\pi RC\sqrt{6}}$$

PROCEDURE:

- Connect the circuit as shown in Fig.
- A 9 V Supply Voltage is given by using Regulated power supply and output is taken from collector of the Transistor.
- Observe the output signal from CRO and note down the output amplitude and time period (T_d).
- Calculate the frequency of oscillations theoretically and verify it practically.
- Repeat the above procedure for different values of R & C.
- Trace the output waveform from CRO.

OBSERVATIONS:

T_d = time period of oscillations

$$f = 1/T_d$$

Table: for determining the frequency of RC phase shift oscillator.

S.No.	Resistance	Frequency (KHz)
1.		
2.		
3.		
4.		
5.		

CALCULATIONS:

The frequency of oscillations of RC-Phase Shift Oscillator is,

$$\frac{1}{2\pi RC\sqrt{6}}$$

RESULT:

C. The frequency of phase shift oscillator is calculated for different values of R and C.

SOURCE OF ERROR AND PRECAUTIONS:

- All connections should be right.
- Readings should be taken properly.
- Power supply should be switched off after use.

MAULANA AZAD
NATIONAL INSTITUTE OF TECHNOLOGY, BHOPAL
DEPARTMENT OF PHYSICS

M.Sc. II Sem

EXPERIMENT: R-C COUPLED AMPLIFIER

AIM:

To study the frequency response of two stage RC coupled amplifier and to determine the bandwidth.

APPARATUS:

Transistor, resistor, capacitor, voltmeter, Power supply, connecting wires.

THEORY:

Description of components: D.C. power supply, the resistances R_1 , R_2 and R_E provides potential divider biasing and stabilization network. i.e. It establishes a proper operating point to get faithful amplification. R_E reduces the variation of collector current with temperature. The potential divider bias provides forward bias to the emitter junction and reverse bias to the collector junction. Since the emitter is grounded, it is common to both input and output signals. Therefore, the amplifier is common-emitter amplifier. Capacitor C_{in} ($= 10 \mu F$) isolates the d.c. component and the internal resistance of the signal generator and couples the a.c. signal voltage to the base of the transistor. The capacitor C_E connected across the emitter resistor R_E is of large value ($= 100 \mu F$) offers a low reactance path to the alternating component of emitter current and thus bypasses resistor R_E at audio frequencies. Consequently, the potential difference across R_E is due to the d.c. component of the current only. The coupling capacitor C_c ($= 10 \mu F$) couples the output of the first stage of amplifier to the input of the second stage. It blocks the d.c. voltage of the first stage from reaching the base of the second stage. The output voltage is measured between the collector and emitter terminals.

When a.c. signal is applied to the base of the first transistor, it is amplified and developed across the out of the 1st stage. This amplified voltage is applied to the base of next stage through the coupling capacitor C_c where it is further amplified and reappears across the out put of the second stage. Thus the successive stages amplify the signal and the overall gain is raised to the desired level. Much higher gains can be obtained by connecting a number of amplifier stages in succession (one after the other). Resistance-capacitance (RC) coupling is most widely used to connect the output of first stage to the input (base) of the second stage and so on. It is the most popular type of coupling because it is cheap and provides a constant amplification over a wide range of frequencies. Fig. 1 shows the circuit arrangement of a two stage RC coupled CE mode transistor amplifier where resistor R is used as a load and the capacitor C is used as a coupling element between the two stages of the amplifier.

$$\text{Voltage Gain} = V_{out}/V_{in}$$

Where

V_o = Output voltage

V_{in} = Input voltage

$$\text{Bandwidth of the amplifier} = f_2 - f_1 \text{ KHz}$$

Where

f_1 = lower cut-off frequency

f_2 = upper cut-off frequency


Fig. 1

PROCEDURE:

- The circuit connections are made as shown in the Fig.1.
- First the signal generator is connected directly to the a.c. milli-voltmeter by keeping signal frequency at about 500 Hz. The amplitude (voltage) of the input signal is adjusted to 0.1V or 0.05V. This is the amplifier input (V_i).
- Now the signal generator is disconnected from the a.c. milli-voltmeter and connected to the input of the amplifier and the a.c. milli-voltmeter is connected to the output of the amplifier.
- Set the input frequency at 10 Hz, note the output voltage (V_o) from the a.c. milli-voltmeter and keeping the input voltage, V_i constant.
- Vary the input frequency ' f ' and note the output voltage.
- The frequency of the input signal is varied in convenient steps i.e. at least 5 values with equal intervals, in each range of frequency in the signal generator, the output voltage V_o is noted in the table for each frequency.
- Calculate the voltage gain, G of the amplifier for each value of the frequency, f of the input signal, using the relation, Voltage gain, $G = V_o / V_{in}$.

OBSERVATIONS:

Maximum frequency

Lower cut-off frequency

Upper cut-off frequency

Input voltage $V_i = \dots\dots\dots V$

Table: for calculate the Voltage gain

S.No.	Applied frequency	Output Voltage	Gain= V_{out}/V_{in}	Gain in dB ($20\log_{10}$ V_{out}/V_{in})
1.				
2.				
3.				
4.				
5.				

To determine the bandwidth (BW) of the amplifier

Draw the frequency response curve as said above, by taking the frequency f (or $\log_{10} f$) on X-axis and voltage gain on Y-axis. Note the maximum gain, G_{max} and mark the value of $0.707G_{max}$ on the y-axis. From that value draw a line (dashed line) parallel to x-axis. This line cuts the curve at two points, called the half-power points. From those two points draw two perpendicular lines on to x - axis, the feet of two perpendiculars corresponding to two frequencies f_1 and f_2 . These are called as lower half power frequency and the upper half-power frequency (or cut-off frequency). The difference between these two frequencies f_1 and f_2 is the bandwidth (BW) of the amplifier.

$\therefore$ Bandwidth of the amplifier = $f_2 - f_1$


Figure 2 Frequency Response Curve of a RC Coupled Amplifier

CALCULATIONS:

$$\text{Voltage Gain} = V_{out}/V_{in}$$

$$\text{Gain in dB} = 20\log_{10} V_{out}/V_{in}$$

$$\text{Bandwidth of the amplifier} = f_2 - f_1 \text{ KHz}$$

Where, f_1 = lower cut-off frequency

f_2 = upper cut-off frequency

RESULT:

The graph is plotted between gain and frequency. This curve is known as frequency response curve and the values of f_2 and f_1 from graph are.....

SOURCE OF ERROR AND PRECAUTIONS:

1. Before going to the experiment the input voltage V_i should be measured.
2. The input voltage should be maintained at constant value throughout the experiment.
3. All connections should be tight.
4. Transistor terminals must be identifying properly.
5. Reading should be taken properly.
6. Scale on the graph should be taken correctly.

MAULANA AZAD
NATIONAL INSTITUTE OF TECHNOLOGY, BHOPAL
DEPARTMENT OF PHYSICS

M.Sc. II Sem

EXPERIMENT: PUSH-PULL AMPLIFIER

AIM:

To study the Output gain, output power and frequency response of Push-Pull amplifier.

APPARATUS:

DC power supply, Two L plate transformers, transistors, resistors and capacitor.

THEORY:

Push-pull amplifier is a power amplifier and is frequently employed in the output stages of electronic circuits. It is used whenever high output power at high efficiency is required. Diagram shows circuit diagram of push-pull amplifier. Two transistors T_{R1} and T_{R2} placed back to back are employed. Both transistor are operated in class B operation i.e. collector current is nearly zero in the absence of the signal. The centre tapped secondary of driver transformer T_{R1} supplies equal and opposite voltage to the base circuit of two transistors. The output transistor T_{R2} has the centre-tapped primary winding. The supply voltage V_{cc} is connected between the bases and this centre tap. The output load is connected across the secondary of this transformer.


Input signal appears across the secondary AB of driver transformer. Suppose during the first half-cycle of the signal, end A becomes positive and end B negative. This will make the base emitter junction of T_{R1} reverse biased and that of T_{R2} forward biased. The circuit will conduct current due to T_{R2} only. Therefore, this half-cycle of the signal is amplified by T_{R2} only and appears in the lower half of the primary of the output transformer. In the next half cycle of the signal, T_{R1} is forward biased whereas T_{R2} is reverse biased. Therefore, T_{R1} conducts consequently this half-cycle of the signal amplified by T_{R1} and appears in the upper half of the output transformer primary. The centre tapped primary of output transformer combines two collector current to form a sine wave output in secondary. It may be noted here that push-pull arrangement also permits a maximum transfer of power to the load through impedance matching.

Voltage gain

$$A_v = V_{out}/V_{in}$$

Output power

$$P = V^2/R \text{ (V is the RMS value of the output signal).}$$

PROCEDURE:

- Connect audio frequency function generator across input sockets and set it at sine wave signal of 20 mV-50mV peak to peak amplitudes, 1KHz frequency.
- Connect CRO across output sockets.
- Switch ON the instrument using ON/OFF toggle switch provided on the front panel.
- Observe the amplified output at CRO. Note down the output amplitude.
- Calculate the voltage gain of amplifier, using formula
$$A_v = V_{out}/V_{in}$$
- Increase the frequency of the signal towards 100 KHz in small steps and note down the voltage gain at different frequencies.
- Note down the observation in table no. (1) and plot a graph between Voltage gain vs Frequency.
- Connect the load resistance (R7 or R8) across output sockets and calculate the output power by using the formula:
$$P = V^2/R$$
 (V is the RMS value of the output signal).

OBSERVATIONS:

- Voltage gain of amplifier, using formula
$$A_v = V_{out}/V_{in}$$
- Output power by using the formula:
$$P = V^2/R$$
 (V is the RMS value of the output signal).

Table: for determine the gain of the signal

S.No.	Frequency	Input Signal	Output Signal	Gain Output/Input
1.				
2.				
3.				
4.				
5.				

CALCULATIONS:

$$A_v = V_{out}/V_{in}$$

$$P = V^2/R$$

RESULT:

The Voltage gain A_v _____

Output Power P = _____

SOURCE OF ERROR AND PRECAUTIONS:

- Connection should be right.
- Readings should be taken properly.
- Supply should be switched-off after use.

Low Pass and High Pass Filters

Aim

To introduce to Low Pass and High Pass Filters. Determine the amplitude and frequency response characteristics of Low and High Pass filters. Draw graphs showing the frequency responses of an low-pass filter (LPF) and an high-pass (HPF) filter.

Apparatus Required

Resistors, Capacitor, CRO, Function Generator, Bread Board, Wires.

Theory

Low-pass and high-pass filters are two commonly used types of filters that work in opposite ways to filter signals. Low-pass filters, as the name suggests, allow low-frequency signals to pass through while attenuating high-frequency signals. On the other hand, high-pass filters allow high-frequency signals to pass through while attenuating low-frequency signals.

Low Pass Filter

The circuit diagram of the low pass filter is shown below. The circuit of LPF can be built with a resistor as well as a capacitor in series so that the output can be achieved. Once the input is given to the circuit of the LPF, then the resistance will give a stable obstacle, however, the capacitor position will have an effect on the output signal.


Figure 1 RC Low Pass Filter


Figure 2 Frequency Response of RC Low Pass Filter

High Pass Filter

The circuits of a high-pass filter and a low-pass filter are completely opposite, because the two components have been swapped, and the filter output signal is now taken out of the resistor

The low-level filter circuit only allows the signal to pass below its cut-off frequency point f_c . As the name suggests, the passive high-pass filter circuit only passes the signal above the selected cut-off point, and eliminates any low-frequency signals in the waveform. Consider the following circuit:


Figure 3 RC high pass filter


Figure 4 Frequency response of high pass filter

Formula:

The cut-off frequency of low pass and high pass filter is given by

$$f_c = 1/2\pi RC$$

Where R = Resistance used in the circuit

C= Capacitance used in the circuit.

Procedure:

1. Construct a circuit like in figure 1 by connecting the capacitor and the resistor in series for low pass filter.
2. Connect the oscilloscope across the capacitor for low pass filter.
3. Calculate the cutoff frequency from $f_c = \omega c/2\pi = 1/2\pi RC$.
4. Draw frequency response curve [$\log_{10} f$ vs $\log_{10} (v_o/v_{in})$] for both low pass and high pass filter.
5. Compare the theoretical cutoff frequency with experimental cutoff frequency.
6. Repeat the procedure for high pass filter.

Input signal Frequency f (Hz)	$\log_{10} f$	V_o (volts)	$\log_{10} (v_o/v_{in})$

Table for frequency response curve for both Low pass and high pass filter

Result:

The frequency response curves of the first-order RC low-pass and high-pass filters were analyzed, and their characteristics were plotted.

MAULANA AZAD NATIONAL INSTITUTE OF TECHNOLOGY, BHOPAL
DEPARTMENT OF PHYSICS

M.Sc. I Sem. / II Sem.

EXPERIMENT: INVERTING & NON INVERTING AMPLIFIER

AIM:

To study gain of Inverting and Non-Inverting amplifier and also study the frequency response characteristics and find out the bandwidth.

APPARATUS:

Operational amplifier as an Inverting and Non-Inverting Amplifier Kit, Patch chords, CRO.

THEORY:

INVERTING AMPLIFIER-

An inverting amplifier using Op-Amp is a type of amplifier using Op-Amp where the output waveform will be phase opposite to the input waveform. The input waveform will be amplified by the factor A_v (voltage gain of the amplifier) in magnitude and its phase will be inverted. In the inverting amplifier circuit the signal to be amplified is applied to the inverting input of the Op-Amp through the input resistance R_i . R_f is the feedback resistor. R_f and R_{in} together determine the gain of the amplifier. Inverting operational amplifier gain can be expressed using the equation $A_v = -R_f/R_i$. Negative sign implies that the output signal is negated. The circuit diagram of a basic inverting amplifier using Op-Amp is shown below.


Figure 1. Inverting Operational Amplifier

The input and output waveforms of an inverting amplifier using Op-Amp is shown below. The graph is drawn assuming that the gain (A_v) of the amplifier is 2 and the input signal is a sine wave. It is clear from the graph that the output is twice in magnitude when compared to the input ($V_{out} = A_v \times V_{in}$) and phase opposite to the input.

A simple practical inverting amplifier using 741 IC is shown in Figure 1. A 741 is a high performance and of course the most popular operational amplifier. It can be used in a variety of applications like Integrator, Differentiator, Voltage follower, Amplifier etc. A

741 has a wide supply voltage range ($\pm 22\text{V DC}$) and has a high open loop gain. The IC has an integrated compensation network for improving stability and has short circuit protection. Signal to be amplified is applied to the inverting pin (pin2) of the IC. Non inverting pin (pin3) is connected to ground. R_i is the input resistor and R_f is the feedback resistor. R_f and R_i together sets the gain of the amplifier.

Mathematically the voltage gain offered by the circuit is given as

$$A_v = \frac{V_o}{V_i}$$

$$\begin{aligned} V_i - V_1 &= I_i R_i \\ V_1 - V_o &= I_f R_f \end{aligned}$$

However, we know that an ideal op amp has infinite input impedance due to which the currents flowing into its input terminals are zero i.e. $I_1 = I_2 = 0$. Thus, $I_i = I_f$. Hence,

$$\begin{aligned} V_i - V_1 &= I_f R_i \\ V_1 - V_o &= I_f R_f \end{aligned}$$

We also know that in an ideal op amp the voltage at inverting and non-inverting inputs are always equal.

As the non-inverting terminal is grounded, zero voltage appears at the non – inverting terminal. That means $V_2 = 0$. Hence, $V_1 = 0$, also. So, we can write

$$\begin{aligned} V_i - 0 &= I_f R_i \\ 0 - V_o &= I_f R_f \end{aligned}$$

from, above two equations, we get,

$$\begin{aligned} -\frac{V_o}{V_i} &= \frac{I_f R_f}{I_f R_i} \Rightarrow \frac{V_o}{V_i} = -\frac{I_f R_f}{I_f R_i} \\ &\Rightarrow \frac{V_o}{V_i} = -\frac{R_f}{R_i} \end{aligned}$$


The voltage gain of the inverting Op-Amp is,

$$A_v = -\frac{R_f}{R_i}$$

This indicates that the voltage gain of the inverting amplifier is decided by the ratio of the feedback resistor to the input resistor with the minus sign indicating the phase-reversal. Inverting amplifiers exhibit excellent linear characteristics which make them ideal as DC amplifiers. Moreover, they are often used to convert input current to the output voltage in the form of Trans-resistance or Trans-impedance Amplifiers. Further, these can also be used in audio mixers when used in the form of Summing Amplifiers.

NON-INVERTING AMPLIFIER

A Non-inverting Op-Amp uses Op-Amp as a main element. The Op-Amp has two input terminals (pins). One is inverting denoted with minus sign (-), and other is non-inverting denoted with a positive sign (+). When we apply any signal to the Non-Inverting input, it does not change its polarity when it gets amplified at the output terminal. So, in that case, the gain of the amplifier is always positive.

Here, in the above circuit, we connect an external resistance R_1 and feedback resistance R_f at inverting input. Now, by applying Kirchhoff Current Law, we get,

$$\frac{v_1}{R_1} = \frac{v_o - v_1}{R_f} \dots \dots \dots (i)$$

Let us assume the input voltage applied to the non-inverting terminal is v_i .

Now, if we assume that the op amp in the circuit is ideal op amp, then,

$$v_1 = v_i$$

Therefore, equation (i) can be rewritten as,

$$\frac{v_i}{R_1} = \frac{v_o - v_i}{R_f}$$

$$\Rightarrow v_i \frac{R_f}{R_1} = v_o - v_i$$

$$\Rightarrow v_o = v_i \left(1 + \frac{R_f}{R_1} \right)$$

$$\Rightarrow \frac{v_o}{v_i} = \left(1 + \frac{R_f}{R_1} \right)$$

The closed loop gain of the circuit is,

$$A = \left(1 + \frac{R_f}{R_1} \right)$$

This term does not contain any negative part. Hence, it proves that the input signal to the circuit gets amplified without changing its polarity at the output.

From the expression of voltage gain of a Non-Inverting Op-Amp, it is clear that, the gain will be unity when $R_f = 0$ or $R_1 \rightarrow \infty$.

When, $R_f = 0$

$$A = \left(1 + \frac{R_f}{R_1} \right) = \left(1 + \frac{0}{R_1} \right) = 1$$

When, $R_1 \rightarrow \infty$

$$A = \left(1 + \frac{R_f}{R_1} \right) = \left(1 + \frac{R_f}{\infty} \right) = 1$$

PROCEDURE

For Inverting/Non Inverting Amplifier

1. Connect the circuit as shown in Figure 1/2.
2. Use R_1 in the input circuit and r_f in the feedback circuit. ($R_1=10K$ and $R_f=10K$)
 $\underline{R_2}$
3. Set the input supply at 1V.
4. Note down the Output using DC Voltmeter.
5. Use R_2 instead of R_1 and note down the output voltage. ($R_2=1K$)
6. Repeat steps (2-4) for different input voltages (0.5V and 1.2 V)

For Frequency Response

1. Connect the circuit as shown in Figure 3 i.e. connect audio frequency sinewave generator across input of the circuit. Set the output of generator to 2V P-P amplitude, 10 Hz frequency.


Figure 3. Frequency Response

2. Connect CRO probes across output of the circuit.
3. Switch ON the instrument using ON/OFF toggle switch provided on the front panel.
4. Gradually increase the frequency of sine wave generator towards 100 KHz and every time note down the output amplitude from CRO. Note down all the observation in Table No. 1. And calculate the gain of amplifier using formula.

$$A_v = \frac{V_o}{V_i}$$

5. Draw a graph between frequency vs gain by taking frequency along X-axis and gain along Y-axis.

Table 1.

S.No.	Frequency	V _{in}	V _{out}	Gain=20log ₁₀ V _{out} /V _{in}

RESULT:

The frequency response of an Op-Amp is studied. The gain of Inverting and Non-Inverting Amplifier is

SOURCE OF ERROR AND PRECAUTIONS:

1. Switch ON the power supply only after making all connections.
2. CRO must be handled carefully.

MAULANA AZAD NATIONAL INSTITUTE OF TECHNOLOGY, BHOPAL
DEPARTMENT OF PHYSICS
M.Sc. I Sem. / II Sem.

EXPERIMENT: INTEGRATOR & DIFFERENTIATOR

AIM: To study frequency response characteristics of Integrator and Differentiator.

APPARATUS:
Operational Amplifier as a Differentiator and Integrator Kit, Patch chords, CRO

THEORY:

INTEGRATOR-

When a voltage, V_{in} is firstly applied to the input of an integrating amplifier, the uncharged capacitor C has very little resistance and acts a bit like a short circuit (voltage follower circuit) giving an overall gain of less than 1, thus resulting in zero output. As the feedback capacitor C begins to charge up, its reactance X_c decreases and the ratio of Z_f/R_1 increases producing an output voltage that continues to increase until the capacitor is fully charged. At this point the ratio of feedback capacitor to input resistor (Z_f/R_1) is infinite resulting in infinite gain and the output of the amplifier goes into saturation. (Saturation is when the output voltage of the amplifier swings heavily to one voltage supply rail or the other with no control in between). The circuit design generate triangular wave providing square wave as input to the integrator. Hence, the integrator circuit generates integral output with respect to the input waveform. The output voltage is directly proportional to the negative integral of the input voltage and inversely proportional to the time constant $R_1 C_F$. For example, if the input is a sine wave, the output will be a triangular wave, as shown in Figure 2. The input offset voltage v_{io} and the part of the input current charging capacitor C_F produces the error voltage at the output of the integrator. Therefore, in the practical integrator shown in Figure 1, a resistor R_F is connected across the feedback capacitor C_F . Thus, R_F limits the low frequency gain and hence minimizes the variations in the output voltage. Both the stability and the low frequency roll off problems can be corrected by the addition of a resistor R_F as shown in Figure 1.


Figure 1. Op-Amp as an Integrator


Figure 2. Input and Output waveforms of an Integrator

DIFFERENTIATOR

A differentiator circuit is one in which the voltage output is directly proportional to the rate of change of the input voltage with respect to time. This means that a fast change to the input voltage signal, the greater the output voltage change in response. As a differentiator circuit has an output that is proportional to the input change, some of the standard waveforms such as sine waves, square waves and triangular waves give very different waveforms at the output of the differentiator circuit.

Figure 3 shows Op-Amp as a differentiator. The differentiator may be constructed from a basic inverting amplifier if an input resistor R_1 is replaced by a capacitor C_1 . Thus the output V_0 is equal to the $R_F C_1$ times the negative instantaneous rate of change of the input V_{in} with the time. Since the differentiator performs the reverse of the integrator's function, a cosine wave input will produce a sine wave output or a triangular input will produce a square wave output. This makes the circuit unstable. Also, the input impedance X_{C_1} decrease with increase in frequency, which makes the circuit very susceptible to high frequency noise. When amplified, this noise can completely override the differentiated output signal. Both the stability and the high frequency noise problems can be corrected by the addition of two components, R_1 and C_F , as shown in Figure 3. This circuit is a practical differentiator.


Figure 3. Op-Amp as a Differentiator


Figure 4. Input and output waveforms of a Differentiator

PROCEDURE

OP-Amp as an Integrator

1. Connect audio frequency function generator of 2V, 1KHz frequency at the input. Select the function generator at the square wave output.
2. Connect the CRO probes at the output sockets.
3. Switch ON both the instruments one by one.
4. Observe the Triangular wave as the output on CRO.
5. Repeat the above steps with different input waves and observe the output on CRO as shown in Figure 2.

OP-Amp as a Differentiator

1. Connect audio frequency function generator of 400 mV and 500 Hz frequency at the input. Select the function generator at the triangular wave output.

2. Connect the CRO probes at the output sockets.

3. Switch ON both the instruments one by one.

4. Observe the Square wave as the output on CRO.

5. Repeat the above steps with different input waves and observe the output on CRO shown in Figure 4.

RESULT:

The Op-Amp as an Integrator and differentiator is studied. The waveforms are traced on graph.

SOURCES OF ERROR AND PRECAUTIONS:

1. Power supply should be switched off after use.
2. Switch the power ON only after making the connections.
3. Handle the CRO carefully.

COMMONLY USED PHRASES

MAULANA AZAD
NATIONAL INSTITUTE OF TECHNOLOGY, BHOPAL
DEPARTMENT OF PHYSICS

M. Sc.. II Sem.

EXPERIMENT: CLIPPING & CLAMPING

AIM:

Study of frequency response characteristics of Clipping & Clamping circuit.

APPARATUS:

Power supply, Clipping & Clamping circuit, Connecting leads, CRO.

THEORY:

Clipping: The clippers have the ability to remove signal voltages above or below a specified level & hence change the wave shape of the I/P signal. Most of the clippers employ diodes & are known as diode clippers.

Different type of clippers are-

Positive & Negative clipper: A circuit that removes +ve half-cycle of the signal is called +ve clipper. Sometimes, it is required to remove the -ve half cycle of the I/P signal, the only thing to be done is to reverse the polarity of the diode connected across load, such a clipper is known as a -ve clipper.

Biased clipper: A clipper used to remove a small portion of +ve or -ve half cycle of the signal Voltage is called a biased clipper. A diode is employed in series with a battery of different volts depending upon the requirement.

Combination clipper: In this circuit small portion of +ve as well as small portion of -ve half- Cycle of the signal voltage is removed.



Figure 1: Clipper circuit

Clamper: A clamping circuit adds d.c component to the signal in such away that it pushes the signal either on the +ve side or on the -ve side. When the circuit pushes the

signal on the +ve side then -ve peak of the signal falls on the zero level, this circuit is called a +ve clamper. When the circuit pushes the signal on the -ve side, this is -ve clamper.



Figure 2: Clamper circuit

PROCEDURE:


CLIPPER:


- Connect the circuit as shown in the circuit diagram
- Observe the output on CRO
- Now connect the battery and observe the output.
- Repeat the above steps for each type of clipper

CLAMPER:

- Connect the circuit as shown in the circuit diagram
- Note down the output on CRO
- Use dc/ac switch of CRO to calculate the added dc level to the wave

WAVEFORM

Positive Clamper

Negative Clamper

RESULT:

Verified the circuit diagram of clipper and clamper waveforms are drawn.

SOURCE OF ERROR AND PRECAUTIONS:

- i. Connection should be proper & tight.
- ii. Switch 'ON' the supply after completing the circuit.
- iii. Handle the equipment with care.
- iv. Loose and wrong connections are to be avoided.
- v. The output waveforms should be obtained without distortion.

M. Sc. II Sem.

EXPERIMENT: LOGIC GATES

AIM:

To study and verify the truth tables of AND, OR and NOT gates.

APPARATUS:

AND gate, NOT gate, OR gate, Power supply, connecting wires.

THEORY:


Logic gates are the digital electronic circuits because they are made up of number of electronic devices and components with one output and one or more inputs. They are the basic building blocks of any logic circuit. AND, OR and NOT gates are basic logic gates. A table which lists all possible combinations of input variables and the corresponding outputs is called a "truth table". It shows how the logic circuit's output responds to various combinations of logic levels at the inputs.

AND GATE:-

An AND gate has two or more inputs but only one output. The logic symbol & truth table of two input AND gate are shown in figure. The symbol for AND operation is ".".

With input variables A & B the Boolean expression for output can be written as

$$Y = A.B$$

INPUT		OUTPUT
A	B	
0	0	0
1	0	0
0	1	0
1	1	1

OR GATE:-

Like an AND gate, an OR gate have two or more inputs but only one output. The logic symbol & truth table of two input OR gate are shown in figure. The symbol for OR operation is "+". With input variables A & B the Boolean expression for output can be written as

$$Y = A + B$$

INPUT		OUTPUT
A	B	
0	0	0
1	0	1
0	1	1
1	1	1

NOT GATE:-

A NOT gate is also known as an inverter, has only one input and only one output. It is a device whose output is always the complement of its input. The logic symbol & truth table of NOT gate are shown in figure. The symbol for NOT operation is "-" (bar). With input variable A the Boolean expression for output can be written as

$$Y = \bar{A}$$

NOT

INPUT	OUTPUT
A	
0	1
1	0

PROCEDURE:

- To verify the truth table of a logic gate, connect inputs and also connect output of logic gate to output indicator.
- Switch ON the instrument using OFF/ON toggle switch provided on the front panel.
- The logical inputs of the truth table are applied and the corresponding output is noted.
- Observe the output indicator. If it glows the indication is that the output is in state '1' and if it does not glow the indication is that the output is in state '0'.
- Similarly the output is noted for all other combinations of inputs as shown in the truth table.
- In this way, the truth table of a logic gate is verified.

RESULT:

The truth table of logic gates AND, OR and NOT is verified.

SOURCE OF ERROR AND PRECAUTIONS:

- All connections should be made neat and tight.
- While making connections main voltage should be kept switched off.
- Never touch live and naked wires.
- Handle the equipment with care.

EXPERIMENT: Rectangular Waveguide

AIM:

To determine the frequency and wavelength in a rectangular waveguide working in TE₁₀ mode.

APPARATUS:

Rectangular Waveguide, Laser Source, Receiver

THEORY:

Any system of conductors and insulators which is used for conveying the electromagnetic wave is called as a wave guide. It is a hollow conducting medium or tube which transmits high frequency electromagnetic waves from source to destination. A specially designed hollow metalloid pipes are used as wave guides. Such constructed wave guides provide an attenuation of transmission line to transmit the electrical energy at microwave frequencies. Inside a wave guide, any configuration of electrical and mechanical fields must have a solution of Maxwell's field equations.

In addition, the boundary conditions imposed by the walls of the guides must have to be satisfied by these fields. There are different modes by which field configurations meet the requirements. These modes include transverse electric or TE modes or H modes. In such modes, electric field is transverse to the axis of the guide and has no component in the direction of the guide except at the location of the associated magnetic field.

In another particular case, magnetic field is everywhere transverse to the guideline axis and some places the electric field has components in the direction of propagation. Such type of mode is called TM or transverse magnetic mode or E mode.

Generally waveguides are classified into two major type's namely Metal waveguides and Dielectric waveguides. These waveguides are in the form of a closed metal pipe. By the reflection from the conducting walls, the wave propagation is characterized inside the metal wave guide. These are used in microwave ovens, microwave radio links, satellite communication, radar sets, etc as a transmission line at microwave frequencies to connect microwave receivers and transmitters to their antennas.


Dielectric waveguides consists of dielectrics only. To propagate the electromagnetic wave along the wave guide, it employs the reflection from dielectric interfaces. Rather than having a hollow pipe, dielectric waveguide employs a solid dielectric rod. In practical, an optical fiber is a dielectric waveguide designed to work at the optical frequencies. Again, these waveguides are classified into two types namely rectangular waveguide and circular waveguide.

In rectangular waveguide, electric field varies with time having maximum at the center whereas the magnetic lines are tangential to the walls.


Both of these are curved round and pass through the guide. The area of this waveguide is in rectangular fashion. These are used at the ultra-high frequencies as alternate to transmission line for transmission of electrical energy.

Circular Waveguide has the uniform circular cross section along their length. These are used as an alternate to the rectangular waveguides. These are often used to feed the conical horns and the round cross section makes it easy to machine. These waveguide supports transverse magnetic (TM) and transverse electric (TE) modes.

Metal Waveguides


Rectangular Waveguide


Circular Waveguide

Dielectric Waveguides


Dielectric Slab Waveguide

Optical Fiber

Figure 1

In transverse electric (TE) modes $E_z = 0$ and $H_z \neq 0$

Fields in the TE mode

The source must operate at a frequency higher than the cutoff frequency to propagate a wave through a waveguide for that particular mode. If the frequency of the waveguide source is less than the cutoff frequency, then the wave is quickly attenuated at the vicinity of the source for that particular mode.

PROCEDURE

1. Switch ON the power supply.
- 2.

RESULT:

The frequency and wavelength in a rectangular waveguide working in TE₁₀ mode are found to be.....

M. Sc. II Sem.

EXPERIMENT: FIELD EFFECT TRANSISTOR (FET)

AIM:

- To Draw the drain and transfer characteristics of a given FET.
- To find the drain resistance (r_d) amplification factor (μ) and Trans conductance (g_m) of the given FET.

APPARATUS:

FET, Regulated power supply, Voltmeter, Ammeter, Bread Board, Connecting Wires.

THEORY:

The Field Effect Transistor is a three terminal unipolar semiconductor device because their action depends only on one type of charge carriers i.e. either electrons or holes. The three terminals are Gate, Drain and source. The current conduction at output is controlled by means of Electric Field (Voltage at the gate) and hence the name is FIELD EFFECT TRANSISTOR. Schematic symbols for P-channel and N-channel FET are shown in Fig.1.


Fig.1 Symbols for P-channel and N-channel FET

Static Characteristics of an FET:

The following two characteristics are

- Drain characteristics:** Drain characteristics are obtained between the drain to source voltage (V_{DS}) and drain current (I_D) taking gate to source voltage (V_{GS}) as the constant parameter.
- Transfer characteristics:** Transfer characteristics are obtained between the gate to source voltage (V_{GS}) and drain current (I_D) taking drain to source voltage (V_{DS}) as the constant parameter.

FET Parameters:

- Drain Resistance (r_d):** It is given by the relation of small change in drain to source voltage (V_{DS}) to the corresponding change in Drain Current (I_D) for a constant gate to source voltage (V_{GS}), when the FET is operating in pinch-off region.

$$r_d = \frac{\Delta V_{DS}}{\Delta I_D} \text{ at constant } V_{GS} \text{ (from drain characteristics)}$$

2. **Trans Conductance (g_m):** Ratio of small change in drain current (I_D) to the corresponding change in gate to source voltage (V_{GS}) for a constant V_{DS} .

$$g_m = \frac{\Delta I_D}{\Delta V_{GS}} \text{ at constant } V_{DS} \text{ (from transfer characteristics)}$$

3. **Amplification factor (μ):** It is given by the ratio of small change in drain to source voltage (V_{DS}) to the corresponding change in gate to source voltage (V_{GS}) for a constant drain current (I_D).

$$\mu = \left(\frac{\Delta V_{DS}}{\Delta I_D} \right) \times \left(\frac{\Delta I_D}{\Delta V_{GS}} \right) = \frac{\Delta V_{DS}}{\Delta V_{GS}}$$

$$\text{i.e. } \mu = r_d \times g_m$$

PROCEDURE:

- 1) All the connections are made as per the circuit diagram.
- 2) To plot the drain characteristics keep V_{GS} constant at 0V.
- 3) Vary the V_{DD} and observe the values of V_{DS} and I_D .
- 4) Repeat the above steps 2, 3 for different values of V_{GS} at -1V and -2V.
- 5) All the readings are tabulated.
- 6) To plot the transfer characteristics, keep V_{DS} constant at 0.5V.
- 7) Vary V_{DS} and observe the values of V_{GS} and I_D .
- 8) Repeat steps 6 and 7 for different values of V_{DS} at 1V and 1.5V.
- 9) The readings are tabulated.
- 10) From drain characteristics, calculate the values of dynamic resistance (r_d) by using the formula $r_d = \frac{\Delta V_{DS}}{\Delta I_D}$
- 11) From transfer characteristics, calculate the value of trans conductance (g_m) by using the formula $g_m = \frac{\Delta I_D}{\Delta V_{GS}}$
- 12) Amplification factor (μ) = dynamic resistance x Trans conductance

$$\mu = \left(\frac{\Delta V_{DS}}{\Delta I_D} \right) \times \left(\frac{\Delta I_D}{\Delta V_{GS}} \right) = \frac{\Delta V_{DS}}{\Delta V_{GS}}$$

$$\text{i.e. } \mu = r_d \times g_m$$

OBSERVATIONS:

DRAIN CHARACTERISTICS:

S. No.	$V_{GS} = 0V$		$V_{GS} = -1V$		$V_{GS} = -2V$	
	$V_{DS}(V)$	$I_D(mA)$	$V_{DS}(V)$	$I_D(mA)$	$V_{DS}(V)$	$I_D(mA)$
1.						
2.						
3.						
4.						
5.						

TRANSFER CHARACTERISTICS:

S. No.	$V_{DS} = 0.5V$		$V_{DS} = 1V$		$V_{DS} = 1.5V$	
	$V_{GS}(V)$	$I_D(mA)$	$V_{GS}(V)$	$I_D(mA)$	$V_{GS}(V)$	$I_D(mA)$
1.						
2.						
3.						
4.						
5.						

MODEL GRAPH:

DRAIN CHARACTERISTICS


Fig. Drain characteristics

TRANSFER CHARACTERISTICS


Fig. Transfer characteristics

CALCULATIONS:

- 1) Dynamic resistance (r_d)

$$r_d = \frac{\Delta V_{DS}}{\Delta I_D}$$

- 2) Trans conductance (g_m)

$$g_m = \frac{\Delta I_D}{\Delta V_{GS}}$$

- 3) Amplification factor (μ) = dynamic resistance x Trans conductance

$$\mu = r_d \times g_m$$

RESULT:

1. The drain and transfer characteristics of a given FET are drawn.
2. The dynamic resistance (r_d), amplification factor (μ) and Trans conductance (g_m) of the given FET are calculated.

SOURCE OF ERROR AND PRECAUTIONS:

- i. Connections must be made according to circuit diagram.
- ii. Make the connection tight.
- iii. Turn off power supply while making connection.
- iv. The three terminals of the FET must be fully identified.
- v. Source and case should be short circuited.
- vi. Voltage exceeding the ratings of the FET should not be applied.

List of Experiments of Physics Laboratory – III

1. X-Ray Diffraction – Determination of lattice parameters of a crystalline solid
2. UV-Vis Spectrophotometer – Determination of absorption coefficient and bandgap
3. Measurement of ac and dc Electrical Conductivity of bulk and thin film
4. Measurement of e/m using Geiger Muller Counter
5. Dielectric Constant and Curie Temperature of Ferroelectric Ceramics
6. Study of Electron Spin Resonance spectrum of paramagnetic substance
7. Study of Nuclear Magnetic Resonance spectrum of magnetic substance
8. Study the Magnetoresistance behaviour of deposited sample
9. Regulated Power Supply using IC-555, IC-723
10. Flip flops a. SR flip-flop b. Clocked SR flip-flop c. JK flip-flop d. Master-slave flip-flop
e. D- flip-flop f. T- flip-flop
11. Design of full adder/ full subtracter.

CHARACTERISTICS OF G.M.COUNTER :-

AIM:-

Geiger Counting System (GCS), GM and a radioactive source.

PRINCIPAL:-

When a gamma ray (or a charged particle) enters the gas filled GM tube, it ionizes the gas inside it and the electric field applied between the electrodes drifts the electrons towards anode. The electrons thus collected at the anode are counted for various applied voltages using Geiger Counting System. A graph is plotted for applied Vs corrected counts ($N - N_b$) and hence the operating voltage is determined from the graph.

PROCEDURE:-

The Geiger Counter System is connected to the GM tube which is mounted on a stand (virtual mount). The radioactive source is placed in the source holder at a distance of about 5cm from the tube. The GCS is switched on and the counter is reset to zero. The high voltage is increased slowly from minimum until the counting just starts. This threshold voltage is noted. The preset time is set to be 20 seconds and the number of counts for this voltage is recorded. Now, the voltage (V) is increased in steps (say 20V) and the number of counts (N) is recorded every time. Increasing the voltage is stopped when the count rate suddenly increases. Any further increase in voltage may damage the GM tube. The number of counts starts decreasing at this point. In this particular case the voltage should not be increased more than 650V.

The voltage checked on digital multimeter at range 1000D.C. after every setting, before taking readings the meter lead is disconnected so as to avoid unnecessary load across the GM tube. By removing the radioactive source the background count (N_b) is recorded for 20 seconds.

A graph is plotted for applied voltage (V) Vs corrected count rate ($N - N_b$). The threshold voltage and limits of the Geiger plateau are marked. The midpoint of the plateau region gives the operating voltage of the tube. The tube must always be operated with this voltage when it is used.

RESULTS:-


(i) Threshold voltage =V

(ii) Operating voltage =V

1. Reset Counter after every set of reading by pressing reset switch.
2. Take out, Radioactive Source after every reading with the help of fork.
3. Do not hold radioactive Source with the fingers directly.

Graph:-


OA = Threshold Voltage
BC = GM Plateau
OD = Operating Voltage

Reading from GCS
Background counter for 20s,

Trial No.	Applied Voltage V Volt	Counts For 20s N	Corrected Counts / 20s (N-N _B)
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
11			
12			
13			
14			
15			
16			
17			
18			
19			

8

INSTRUCTION FOR HANDLING END WINDOW:-

1. Hold the detector with body in the middle only (when used bare).
2. Do not touch the end window with finger. You may rupture the end window thus damaging it permanently.
3. Carefully take out the detector from the packing box and GM tube Stand with help of detector holder given with the stand.
4. Now connect the detector Red Socket with Red Socket of GM COUNTER and Black with Black.
5. Do not interchange or reverse connect these leads. If you reverse connect these leads you will damage the detector permanently.
6. Once mounted into the stand, do not remove frequently. After use, keep the G.M stand with detector together in a safe place.
7. Do not operate the detector in the DISCHARGE for a long time. Life will be reduced and detector may fail if you operate for long time.
8. **One year warranty against manufacturing defect, Breakage in tube & not working (Due to negligence/weather condition are not covered under warranty.)**

**RADIOACTIVE
SOURCE**

INSTRUCTIONS FOR HANDLING RADIOACTIVE SOURCE:-

1. Carefully take out the source from the packing box. Keep a safe distance from it.
2. Hold the source with body (Container's). Do not touch it with finger. Wash your hand with soap if you do so.
3. Activity, Dose rate and date of manufacturing is mentioned on the back side of source.
4. The source supplied with GM COUNTER must not be resold or transferred to any other user or disposed of except as advised by AERB/RSO.
5. The institution has to provide facilities for safe/Locked Storage, handling and use of radioactive material prescribed by the department of Atomic Energy Regulatory.
6. Source and Lead Containers, must be returned to the manufacturer when not in use. Do not sell any metal parts supplied with instrument to any scrap vendor/kbariwalla in any conditions.
7. Any complaint or discrepancy observed should be brought to your RSO notice within 2 working days.

Aim: To determine the dielectric constant of solid and liquid.

Apparatus Used: Dielectric Constant of Solids & Liquids, Sample Jar (Glass), Sample (Solid): Glass Plate, Bakelite Sheet, Teflon, PZT (Lead Zirconate Titanate), Sample (Liquid): Carbon Tetrachloride (500ml).

Theory:

Dielectric or electrical insulating materials are the substances in which electrostatic field can persist for long times. When a dielectric is placed between the plates of a capacitor and the capacitor is charged, the electric field between the plates polarizes the molecules of the dielectric (Figure 1). This produces concentration of charge on its surface that creates an electric field which is antiparallel to the original field (which has polarized the dielectric). This reduces the electric potential difference between the plates. Considered in reverse, this means that, with a dielectric between the plates of a capacitor, it can hold a larger charge. The extent of this effect depends on the dipole polarizability of molecules of the dielectric, which in turn determines the dielectric constant of the material.


Figure 1

The electrons in the molecules shift toward the positively charged plate. The molecules then create a leftward electric field that partially annuls the field created by the charged plates. (The air gap is shown for clarity; in a real capacitor, the dielectric is in direct contact with the plates.)

The method for determination of dielectric constants of liquids consists in the successive measurement of capacitance, first in a vacuum, and then when the capacitor is immersed in the liquid under investigation. A cylindrical capacitor has been used for liquid samples and different size parallel plate capacitors for solid samples.

Parallel Plate Capacitor:

The capacity of a parallel plate capacitor is given by

$$C = k\epsilon_0 A/d \quad \text{farad,} \quad (1)$$

$$\text{farad,} \quad (1)$$

where A is plate area and d is distance between the plates. It is assumed that the dielectric completely fills the space between the capacitor plates.

Cylindrical Capacitor:

The capacitance per unit length of a long cylindrical capacitor immersed in a medium of dielectric constant k is given by

$$C' = \kappa \frac{2\pi\epsilon_0}{\ln(r_2/r_1)} \text{ farad/m.} \quad (2)$$

Here $\epsilon_0 = 8.854 \times 10^{-12} \text{ C}^2\text{m}^{-2}\text{N}^{-1}$ is permittivity of free space, r_1 is external radius of the inner cylinder and r_2 is internal radius of the outer cylinder.

In actual practice, there are errors due to stray capacitances at the ends of the cylinders and the leads. In any accurate measurement, it is necessary to eliminate these. It has been done in the following way:


Figure 2a: Coaxial Cylinder Set-up (CCS)


Figure 2b
Dielectric Measurement in Solids & Liquids (Non Conducting)

Consider a cylindrical capacitor of length L (in meter) filled to a height $h < L$ with a liquid of dielectric constant κ . Its total capacitance is given by

$$C = \frac{2\pi\epsilon_0}{\ln(r_2/r_1)} [\kappa h + \kappa_0 (L - h)] + C_s \quad (3)$$

Here κ_0 is the dielectric constant of air and C_s is sum total of stray capacitances. A simplification of Eq.(2) leads to

$$C = \frac{2\pi\epsilon_0}{\ln(r_2/r_1)} (\kappa - \kappa_0) h + C_0, \quad (4)$$

where


$$C_0 = C_s + \frac{2\pi\epsilon_0}{\ln(r_2/r_1)} \kappa_0 L.$$

Eq. (3) shows that the measured capacity C is a linear function of h (the height upto which the liquid is filled in the capacitor). If we vary the liquid height h , and measure it, together with the corresponding capacitance C , the plot of the data should be a straight line. The slope

$$\frac{dC}{dh} = \frac{2\pi\epsilon_0}{\ln(r_2/r_1)} (\kappa - \kappa_0) \text{ farad/m} \quad (5)$$

will yield a value for κ if κ_0 , r_1 and r_2 are known. The uncertainty due to C_s has thus been eliminated.

Dilelectric Constant of Solid & Liquids, Model DSL-01
Sample: Carbon Tetrachloride (CCl_4)

EXPERIMENTAL PROCEDURE

SETUP:

Choose the appropriate Sample Holder, as per sample size and connect the same to the Digital Capacitance Meter. Now add a little space among two plates of the sample holder and switch the capacitance meter 'ON' without keeping any sample. Press "ZERO ADJ" switch for 2-3 seconds, until the meter shows "Zero Cal". The capacitance meter is calibrated now.

SOLID SAMPLE:

The upper plate of the appropriate parallel plate capacitor of the arrangement for solids is raised by turning the top screw anticlockwise to keep/ insert the sample. The plate is then lowered by turning the top screw clockwise till the capacitor plates touch the sample surface. Switch the Meter "ON" and record reading.

LIQUID SAMPLE:

The glass container is put on the stand below the Coaxial Cylinder Set-up (CCS) as shown in Figure 2b. The container is filled with the experimental liquid up to a height of about 12 cm. The CCS is now lowered until the liquid touches the zero mark on the scale. The Capacity Meter is now adjusted to zero value by pushing zero set button for about a second (pressing it for longer duration will automatically set it to Auto Calibration mode, which is not desirable for this experiment). Remember, the result depends not on the actual value of capacitance, but on its variation with the height of the liquid. The CCS is then lowered in steps of 1 cm and the readings on the Capacity Meter are taken and recorded.

TYPICAL OBSERVATIONS TABLE, CALCULATIONS

SOLID SAMPLE:

- (1) Room temperature: 26°C
- (2) Diameter of smaller capacitor plates: 10.0 mm
- (3) Diameter of bigger capacitor plates: 50.0 mm
- (4) Thickness of lead zirconate titanate Pellet: 2.55 mm
- (5) Thickness of bakelite sheet sample = 1.5 mm
- (6) Thickness of glass plate sample = 1.5 mm
- (7) Thickness of teflon sample = 1.0 mm

Capacity Measurement

Sr. No.	Sample	C (pF)
1	PZT (Lead Zirconate Titanate)	
2	Glass	
3	Bakelite	
4	Teflon	

Liquid Sample:

Sample: Carbon Tetrachloride (CCl_4)

Dielectric Constant of free air (ϵ_0): 1.0059

External radius of the inner cylinder: 25.4 mm

Internal radius of the outer cylinder: 30.6 mm

Capacity Measurement:

S.No.	Liquid Level (cm)	Capacity (pf)
1	0.0	
2	1.0	
3	2.0	
4	3.0	
5	4.0	
6	5.0	
7	6.0	
8	7.0	
9	8.0	

1. A plot of Capacitance v/s Liquid height is shown in Fig. 3.
2. The graph is linear. Its slope dC/dh is measured.

Using Eq.(4),

$$\kappa - \kappa_0 = \frac{(dC/dh)}{2\pi\epsilon_0} \ln(r_2/r_1)$$

One can calculate dielectric constant.

PRECAUTIONS AND SOURCES OF ERROR

1. The connecting leads internal capacitance through compensated by Zero Setting, but still could introduce some error.
2. Presence of vapour over the surface of the liquid.

Aim: To study the biological free radicals for quantitative and qualitative analyses of reactive oxygen species (ROS) and reactive nitrogen species (RNS)

Apparatus used: ESR set up, Helmholtz coil, oscillator

Theory:

1. **Basic Circuit:** The first stage of the ESR circuit consists of a critically adjusted (marginal) radio frequency oscillator having a frequency range of approximately 12–16MHz. A marginal oscillator is required here so that the slightest increase in its load decreases the amplitude of oscillation to an appreciable extent. The sample is kept inside the tank coil of Helmholtz coils. At resonance, i.e. when the frequency of oscillation equal to the Larmor's frequency of the sample, the oscillator amplitude registers a dip due to the absorption of power by the sample. This obviously, occurs periodically - four times in each complete cycle of the Helmholtz coils supply voltage. The result is in amplitude modulated carrier (Fig. 4A) which is then detected using a diode detector and amplified by a chain of three low noise, high gain audio - frequency amplifiers of excellent stability. A sensitivity control is provided in the amplifier to suit the input requirement of any oscilloscope. A block diagram of the ESR Spectrometer is given in Fig. 1


Figure 1. Block Diagram of the ESR Set up

2. **Phase Shifter:** In order to make it possible to use an ordinary displaying type oscilloscope, instead of a measuring oscilloscope which preserve the phase between X and Y plates signals, a phase shifter is provided. This can compensate the phase difference which is introduced in the amplification stage of the ordinary oscilloscope.

The circuit diagram of the phase shifter is shown in Fig.4(B). The primary of the transformer is fed from the 220V, 50Hz (or 110V, 60Hz) mains and the secondary is centre tapped developing V_1-0-V_1 (say). The operation of the circuit may be explained with the help of the vector diagram shown in Fig.4(B). The vectors OA and BO represent the voltage developed in the secondary, in phase and magnitude. The current flowing in the circuit ADB leads the voltage vector BA due to the presence of capacitor C and is shown in the diagram as

I. Voltage developed across resistance R , i.e. V_R is in phase with the current I , and the voltage across capacitor V_C is 90° (lag) out of phase with the current. The vector sum of V_C and V_R is equal to $2V_1$. These are also plotted in the diagram. It is clear from the diagram that as R is varied, V_R will change and the point D will trace a semicircle, shown dotted. The vector OD , or the voltage across points O and D , will, therefore, have a constant magnitude equal to V_1 and its phase, variable from 0 to 180° . This is the voltage which is fed to the X-amplifier of the oscilloscope to correct for any phase change which might have taken place in the rest of the circuit.

3. 50 Hz Sweep Unit: For modulation with a low frequency magnetic field, a 50 Hz current flows through the helmholtz coils. As the resonance in this frequency range occurs at low magnetic fields, no static D.C. magnetic field is required.

4. Power Supplies:

a) D.C. Power Supply: The ESR circuit requires a highly stabilised almost ripple free voltage. These are obtained using integrated circuit regulator.

b) Helmholtz Coils Power Supply: The Helmholtz coils power supply consists of a step-down transformer (220 to 35 V AC). Variable coil current is provided in 10 steps using a band switch, while the current is displayed on a 3 $\frac{1}{2}$ digit panel meter. The output is taken from the two terminals provided on the panel.

5. Helmholtz Coils: There are two coils exactly alike and parallel to each other, so connected that current pass through them in the same direction. The two coils increase the uniformity of the field near the centre.

Number of turns	:	500 in each coil
Diameter of the Windings	:	15.4 cm
Separation of the coil	:	7.7 cm

In the centre of the coils, an attachment is provided to keep the sample in place and to minimise shocks and vibrations.

6. Test Sample: A test sample, Diphenyl Picryl Hydrazyl (DPPH) (Fig. 5) is placed in a plastic tube, which itself is in the induction coils. This increases the filling factor to the maximum. DPPH is a free radical and widely used as a standard for ESR measurements.

8. Controls & Terminals:

- (1) Mains : To switch 'ON' or 'OFF' the ESR Spectrometer.
- (2) Phase: To adjust the phase between X and Y plates signals.
- (3) Current : To control current in Helmholtz coils.
- (4) 'H' Coils : Terminals and switch for Helmholtz coils.
- (5) Frequency : To adjust the frequency of the Oscillator.
- (6) X,Y,E : For X, Y and Earth terminals of the Oscilloscope.

9. Oscilloscope: As the Oscilloscope is not supplied with the spectrometer, it is presumed that the worker is already familiar with the control knobs and functioning of the Oscilloscope in question. Any Oscilloscope, normally available in the laboratory of the following specifications or better, will be quite suitable for the observation of ESR resonance :

Screen diameter : 12.5 cm
Vertical amplifier sensitivity : 50 mV/cm

Procedure:

- (1) Switch on 'H. COIL' power and adjust the current at 150 mA.
- (2) Set the front panel controls of ESR Spectrometer as follows
Frequency : Centred
Phase : Centred
- (3) Observe four peaks on the Screen of CRO. Now adjust the FREQUENCY of the Spectrometer and SENSITIVITY of the CRO to obtain the best results (i.e. sharp peaks and good signal to noise ratio).
- (4) Adjust the PHASE knob to coincide the two peaks with the other two as far as possible.
- (5) Adjust the orientation of Helmholtz coils with respect to the main unit for best overlap of base lines.

Calculation:

From equation

$$h\nu_1 = g \mu_0 H$$

or
$$g = \frac{h \nu_1}{H_0 \mu_0}$$

Substitute the measured values of H_0 and ν_1 universal constants $h=6.625 \times 10^{-27}$ erg.sec and $\mu_0=0.927 \times 10^{-20}$ erg/gauss to get the g factor.

Sources of errors and precautions:

1. Experiment should be set up at a quiet place free from mechanical and electrical disturbances.
2. Y - output from the ESR Spectrometer should be taken through shielded cable to minimise external pick-ups.
3. X and Y plates sensitivities of the Oscilloscope should be adjusted such that they should be in the linear range only.
4. Since the Helmholtz current is unstabilised, care should be taken that it should remain constant during the observation.
5. High currents (~ 200 mA) should not be allowed to flow through the Helmholtz coils for an extended period. This will unnecessarily heat the coils which may get damaged after some time.
6. If the peaks do not coincide on the x-scale, check the sinusoidal wave form of the mains voltage which may be distorted due to overloading of the main line due to other heavy gadgets working on the same line.

Aim: To determine the structure of organic molecules in solution and study molecular physics and crystals

Apparatus Used: NMR set up, solutions

Theory:

Nuclear Magnetic Resonance (NMR) is a powerful analytical tool. It was discovered by Bloch and Purcell in 1945. NMR spectroscopy is nowadays an important standard method in chemistry and biology. In medical applications it is known as magnetic resonance imaging (MRI) and supplements other scanning techniques.

It is based on the magnetic properties of the atomic nucleus. A nucleus which possesses a nuclear spin I has a magnetic moment μ given by

$$\mu = g \mu_n I \dots\dots\dots (1)$$

Here μ_n is nuclear magneton and g is the nuclear g -factor. The nuclear spin can have $2I + 1$ possible orientation with respect to a magnetic field. In the absence of magnetic field, these orientations are of equal energy. If a magnetic field B is applied or is present at the site of the nucleus, then the energy levels split. This is analogous to the Zeeman Effect. Each level is characterized by a magnetic quantum m and has energy given by

$$\begin{aligned} E_m &= - \mu \cdot B \\ &= - g \mu_n I \cdot B \\ &= - g \mu_n B m, \dots\dots\dots (2) \end{aligned}$$

with m taking values $-I, -(I-1), -(I-2), \dots\dots\dots (I-1), I$. A nucleus with spin $\frac{1}{2}$ will have two possible orientations, Fig.1.

Consider the specimen having these nuclei. When they are in a magnetic field, the initial populations of the energy levels are determined by thermodynamics, as described by the Boltzmann distribution. This is very important, and it means that the lower energy level will contain slightly more nuclei than the higher level. The actual differential between the levels depends on μ , B and the temperature T . It is possible to excite these nuclei into the higher level with electromagnetic radiation. The frequency of radiation needed is determined by the difference in energy between the energy levels. This energy difference is given by

$$\Delta E = g \mu_n B, \dots\dots\dots (3)$$

as $\Delta m = \pm 1$.

In order to understand the absorption of radiation by a nucleus (of spin $\frac{1}{2}$) in a magnetic field, imagine that it is in the lower energy level. Its magnetic moment is aligned with the field and is not opposed to it. Remember that the nucleus also has an angular momentum (the spin) and a torque τ , which is perpendicular to the angular momentum I , acts on it,


Fig. 1


Fig. 2

$$\begin{aligned}\vec{r} &= \mu \vec{r} \times B \\ &= g \mu_n I \times B\end{aligned}\quad \dots (4)$$

In the presence of a magnetic field, the axis of rotation will *precess* around the magnetic field, Fig 2. The frequency of precession is called the Larmor frequency. It is given by

$$\nu = \frac{g \mu_n}{h} B \quad \dots (5)$$

It is identical to the transition frequency $\Delta E/h$. A resonance absorption of electromagnetic power occurs when photons (of the radio-frequency field) bombarding the specimen have the proper energy to excite transitions between these levels, i.e., satisfy the condition of resonance. The strength of the absorption depends upon the difference in population of the levels involved.

When the energy is absorbed by the nucleus, then the angle of precession will change. For a nucleus of spin $1/2$, the absorption of radiation "flips" the magnetic moment so that it now opposes the field, Fig 3.

It is important to realize that only a small proportion of "target" nuclei are in the lower state (and can absorb radiation).

There is now the possibility that by exciting these nuclei, the populations of the higher and lower energy levels will become equal. When this occurs, then there will be no further absorption of radiation. The spin system gets saturated. In order to have a sustained absorption, there must be relaxation processes which return nuclei to the lower energy state.

There are two major relaxation processes: (1) Spin-lattice relaxation, (2) Spin-spin relaxation. The later one is not very effective. Consider the former. The atoms containing these nuclei are in rotational and vibration motion. The magnetic fields generated by these motions interact with nuclear precession and cause the nuclei to lose energy and return to the lower state. The energy that a nucleus loses increases the amount of vibration and rotation within the lattice resulting in a little increase in the temperature of the sample.

Procedure:

The sample under investigation is placed in an induction coil, which is the component of the tank circuit of the oscillator (generator). This is the Zavoisky's technique. It is based on the fact that under certain conditions such as absorption of power from generator, the watt loads (Δw) on the generator changes. This change of Δw is proportional to the change in base current ΔI_b or collector current ΔI_c of the generator. The proportionality, however, holds only the power dissipated by the sample due to absorption is small in comparison to the total losses in the circuit. This change in ΔI_c is detected with the conventional circuits. To make the detection simple and more sensitive, the magnetic field and hence the Larmor frequency of the sample is modulated with a low frequency field 50 Hz in the present set-up.


Fig. 3 Block Diagram of NMR Spectrometer

Procedure

- Switch 'ON' power switch of Unit - 3.
- Switch 'ON' switch of 'H-COILS' and adjust current to 80 mA.
- Insert the sample tube containing preferably Cupric Sulphate from the top of Unit - 2.
- Switch 'ON' Battery switch of Unit - 1.
- Adjust the 'Y input' amplitude knob to 50 mV/div. for proton and 20 mV/div for fluorine.
- Rotate frequency knob (Unit - 1) slowly until the peaks are seen.
- Further slowly rotate this knob till the distance between peaks in 10 ms.
- Note the value of frequency displayed on Unit - 3.
- Please note that the magnetic field at sample coil is 5.02 KG or 0.502 Tesla.

Observation:

- Frequency displayed on unit when the distance between peaks 10 ms for solution at 20.82 MHz
- Frequency displayed on unit when the distance between peaks 10ms for solution: 19.59 MHz

Calculation

From Eq

$$g = \frac{h}{\mu_n} \frac{\nu}{B_o}$$

where $h = 6.625 \times 10^{-34} \text{ Js}$

$\mu_n = 5.051 \times 10^{-27} \text{ J/T}$

User's Manual

MAGNETORESISTANCE SET-UP

Model: DMR-01

Manufactured by:

SCIENTIFIC EQUIPMENT & SERVICES

358/1, New Adarsh Nagar, Roorkee-247 667

Ph.: 01332-272852, Fax: 277118

Email: ses@sestechno.com

Website: www.sestechno.com

Measurement of Magnetoresistance of Semiconductors

It is noticed that the resistance of the sample changes when the magnetic field is turned on. The phenomenon, called magnetoresistance, is due to the fact that the drift velocity of all carriers is not same. With the magnetic field on; the Hall voltage $V = E_y t = |\mathbf{v} \times \mathbf{H}|$ compensates exactly the Lorentz force for carriers with the average velocity; slower carriers will be over compensated and faster one undercompensated, resulting in trajectories that are not along the applied field. This results in an effective decrease of the mean free path and hence an increase in resistivity.


Here the above referred symbols are defines as: $\mathbf{v}$ = drift velocity; E = applied electric field; t =

thickness of the crystal; H = Magnetic field

Experimented Set-up for Magnetoresistance

The set-up consists of the following:

1. Four probe arrangement
2. Sample: (Ge: p-type)
3. Magnetoresistance set-up, DMR-11
4. Electromagnet, EMU-50V
5. Constant Current Power Supply, DPS-50
6. Digital Gaussmeter, DGM-102

(1) Four Probe arrangement

It consists of 4 collinear, equally spaced (2mm) and individually spring loaded probes mounted on a PCB strip. Two outer probes for supplying the constant current to the sample and two inner probes for measuring the voltage developed across these probes. This eliminate the error due to contact resistance which is particularly serious in semiconductors. A platform is also provided for placing the sample and mounting the Four Probes on it.

(2) Sample

Ge Crystal (n-type) dimensions : 10 x 10 x 0.5mm.

(3) Magnetoresistance Set-up, Model DMR-11

This unit consists of a digital millivoltmeter and constant current power supply. The voltage and probe current can be read on the same digital panel meter through a selector switch.

(a) Digital Millivoltmeter

Intersil 3½ digit single chip ICL 7107 have been used. Since the use of internal reference causes the degradation in performance due to internal heating an external reference have been used. Digital voltmeter is much more convenient to use, because the input voltage of either polarity can be measured.

Specifications

Range : 0-200mV (100μV minimum)
Accuracy : ±0.1% of reading ± 1 digit


(b) Constant Current Power Supply

This power supply, specially designed for Hall Probe, provides 100% protection against crystal burn-out due to excessive current. The supply is a highly regulated and practically ripple free dc source.

Specifications

Current : 0-20mA Resolution : 10μA
Accuracy : ±0.2% of the reading ±1 digit
Load regulated : 0.03% for 0 to full load
Line regulation : 0.05% for 10% variation

SCHEMATIC DIAGRAM OF DIGITAL PANEL METER

TEST RESULT OF MAGNETORESISTANCE PROBE S.No. 151

Observation and Calculations

(I) Calibration of EMU-50 at a suitable air-gap ≈ 20 mm

S.No.	Current (A)	Magnetic Field (KG)
1	Min	0.087
2	0.25	0.310
3	0.50	0.537
4	0.75	0.793
5	1.00	1.035
6	1.25	1.290
7	1.50	1.544
8	2.00	2.050
9	2.50	2.550
10	3.00	3.050
11	3.50	3.530
12	4.00	4.000

(II) Magnetoresistance of a Ge-crystal (n-type)

Probe Current $I = 4.00$ mA (Constant for the whole set of readings)

$$\frac{\Delta R}{R} = \frac{R_m - R}{R}, \text{ where Sample Resistance (R) without magnetic Field} = 43.92$$

SNo	Current (A)	Mag. Field H (KG)	Voltage V_m (mV)	$R_m = \frac{V_m}{I}$ (Ω)	$\frac{\Delta R}{R} \times 10^{-3}$	Log ($H \times 10^{-2}$) (KG)	Log $\left(\frac{\Delta R}{R} \times 10^{-3} \right)$
1	Min	0.087	175.7	43.925	0.00	0.94	∞
2	0.25	0.310	175.7	43.925	0.00	1.49	∞
3	0.50	0.537	175.8	43.950	0.57	1.73	-0.24
4	0.75	0.793	175.9	43.975	1.14	1.90	0.06
5	1.00	1.035	176.1	44.025	2.28	2.01	0.36
6	1.25	1.290	176.3	44.075	3.41	2.11	0.53
7	1.50	1.544	176.6	44.150	5.12	2.19	0.71
8	2.00	2.050	177.2	44.300	8.54	2.31	0.93
9	2.50	2.550	178.0	44.500	13.09	2.41	1.12
10	3.00	3.050	178.9	44.725	18.21	2.48	1.26
11	3.50	3.530	179.9	44.975	23.90	2.55	1.38
12	4.00	4.000	180.9	45.225	29.60	2.60	1.47


Nature of Graph: H Vs. $\frac{\Delta R}{R}$ as per sheet attached

Q.C. Engineer: Sunil Sharma


Dated : 10/12/2008

Doc. No. QF/QC/TR/DMR-01/02 ; Revision : 00

MAGNETORESISTANCE OF A Ge CRYSTAL #151 **$\Delta R/R$ vs. Magnetic Field (H)**

MAGNETORESISTANCE OF A Ge CRYSTAL #151 **Log-Log Plot $\Delta R/R$ vs. Magnetic Field (H)**

EXPERIMENT NO: 5**DATE****IC555 TIMER – MONOSTABLE & ASTABLE MULTIVIBRATOR CIRCUITS****A) MONOSTABLE MULTIVIBRATOR**

AIM: To design a Monostable Multivibrator using IC555 and compare it's theoretical and practical pulse width.


APPARATUS: Bread Board.

CRO
Probes
Connecting wires
555 Timer, Resistors, Capacitors


THEORY:

Monostable multivibrator is also called as one-shot Multivibrator. When the output is low, the circuit is in stable state, transistor T1 is ON and Capacitor C is shorted to the ground. However, upon application of a negative trigger pulse to Pin-2, transistor T1 is turned OFF, which releases short circuit across the external capacitor and drives the output High. The capacitor C now starts charging up toward VCC through R. However when the voltage across the external capacitor equals $\frac{2}{3} V_{CC}$, upper comparator's output switches from low to high which in turn derives the output to its low state. And the output of the flip flop turns transistor T1 ON, and hence the capacitor C rapidly discharges through the transistor. The output of the Monostable remains low until a trigger pulse is again applied. Then the cycle repeats. The time during which the output remains high is given by

$$t_p = 1.1 R C$$


Waveforms for IC555 Monostable Multivibrator



Once triggered, the circuit's output will remain in the high state until the set time t_p elapses. The output will not change its state even if an input trigger is applied again during this time interval t_p .

DESIGN:

1. Choose a desired pulse width, say $t_p = 1.1$ ms.
2. Choose a value for capacitor C ($0.1 \mu\text{F}$) and then calculate the value of R by using the equation for t_p .

CIRCUIT DIAGRAM:

PROCEDURE:

1. Connect the components/equipment as shown in the circuit diagram.
2. Switch ON the power supply.
3. Connect function generator at the trigger input.
4. Connect channel-1 of CRO to the trigger input and channel-2 of CRO to the output (Pin 3).
5. Using Function Generator, apply 1 KHz square wave with amplitude of approx. equal to 9 V_{pp} at the trigger input.
6. Observe the output voltage with respect to input and note down the pulse width and amplitude.
7. Now connect channel-2 of CRO across capacitor and observe the voltage across the capacitor and note it down.
8. Compare the practical pulse width noted in the step above with its theoretical value ($t_p = 1.1 RC$)

CALCULATIONS:


THEORETICAL Pulse width

R = C =

$t_p = 1.1 RC =$

PRACTICAL Pulse width

$t_p =$

EXPECTED WAVEFORMS:**RESULT:****QUESTIONS:**

1. What is the other name for monostable multivibrator (MSMV)?
2. When MSMV is in stable state, what is the output level?
3. Why trigger is required in the case of MSMV?
4. Which type of trigger pulse is required for MSMV?
5. What is the formula for the output pulse width of MSMV?
6. How long MSMV stays in unstable state?

(B) ASTABLE MULTIVIBRATOR

AIM: To design an Astable Multivibrator using IC555 and compare it's theoretical and practical time period and duty cycle.

APPARATUS: Bread Board.

CRO

Probes

Connecting wires

555 Timer, Resistors, Capacitors

THEORY:

An Astable multivibrator, often called a free-running Multivibrator, is a rectangular-wave-generating circuit. Unlike the Monostable multivibrator, this circuit does not require an external trigger to change the state of the output, hence the name free running. However, the time during which the output is either high or low is determinate by the Two resistors and a capacitor, which are externally connected to the 555 timer.

Figure 1 shows the 555 timer connected as an Astable multivibrator. Initially, when the output is high, capacitor C starts charging towards V_{cc} through R_A and R_B . However as soon as voltage across the capacitor equals $2/3 V_{cc}$, comparator 1 triggers the flip-flop, and the output switches low. Now the capacitor C starts discharging through R_B and the transistor Q 1. When the voltage across C equals $1/3 V_{cc}$, comparator 2's output triggers the flip-flop, and the output goes high. Then the cycle repeats. The output voltage and the capacitor voltage waveforms are shown in the following figures.

As shown in this figure, the capacitor is periodically charged and discharged between $2/3 V_{cc}$ and $1/3 V_{cc}$, respectively. The time during which the capacitor charges from $1/3 V_{cc}$ to $2/3 V_{cc}$ is equal to the time the output is high and is given by

$$t_c = 0.69 (R_A + R_B) C \quad (1)$$

Similarly, the time during which the capacitor discharges from $2/3 V_{cc}$ to $1/3 V_{cc}$ is equal to the time the output is low and is given by



$$t_d = 0.69 (R_B)C \quad (2)$$

Thus the total time period of the waveform is

$$T = t_c + t_d = 0.69(R_A + 2R_B) \quad (3)$$

Therefore the frequency of oscillation is $f_o = 1/T = 1.45/(R_A + 2R_B)C$

And $\% \text{ Duty cycle} = (t_c/T) * 100 \quad (4)$

CIRCUIT DIAGRAM:**EXPECTED WAVE FORMS:****PROCEDURE:**

1. Connect the components/equipment as shown in the circuit diagram.
2. Switch ON the power supply.
3. Connect channel-1 of CRO to the output (Pin 3).
4. Observe the output voltage and note down the time period and duty cycle.
5. Now connect channel-2 of CRO across capacitor and observe the voltage across the capacitor and note it down.
6. Compare the practical time period and duty cycle.

CALCULATIONS:

THEORETICAL time periods

$$t_c = 0.69 (R_A + R_B) C$$

$$t_d = 0.69 (R_B)C$$

Total time period of the waveform, $T = t_c + t_d$

$$\% \text{ Duty Cycle} = (t_c / T) * 100$$

PRACTICAL (from output waveforms)

time period, $T =$

% Duty cycle =

RESULT:

\

QUESTIONS:

1. What is the other name for Astable multivibrator (AMV)?
2. What is the formula for the time period of the waveform of AMV?
3. What is the formula for the % of Duty cycle?

INSTRUCTION MANUAL

FOR

STUDY OF RS, D FLIP FLOP USING NAND GATE & JK, D FLIP FLOPS USING TTL IC's

Flip Flop circuits using NAND Gate and TTL IC's has been designed to verify the truth table of following :

1. RS & D type Flip Flops using NAND Gates.
2. D type Flip Flop using TTL IC.
3. JK Flip Flop using TTL IC.

The instrument comprises of the following built in parts :

1. Fixed output DC regulated Power supply of 5V.
2. 1Hz Monoshot Clock pulse with pulser switch is provided on the front panel.
3. Four Logic inputs logic '0' & logic '1' selectable using SPDT switches are provided on the front panel.
4. Two red output indicators are also provided on the front panel.
5. IC 7400, 7410, 7474 & 7476 are mounted on the front panel & important connections are brought out on sockets.

THEORY

A flip flop is an electronic circuit that has two stable states, one representing a binary '1' and the other binary '0'. If put into one state, the flip flop will remain in that state as long as power is applied or until it is changed. It thus remembers the data. In digital circuits, flip flops are used in a variety of storage, counting, sequencing and timing applications. There are three basic types of flip flops, the set - reset (also known as 'R-S flip flop or a latch), the 'D' type and the 'JK'. The 'RS' flip flop is the simplest. It has two inputs 'S' & 'R' and two outputs 'Q' and 'Q' . Applying appropriate logic signals to either 'S' or 'R' input will put the latch into one state

or the other. When a flip flop is set by 'S' input, it is said to be storing binary '1' ('Q' output is = High). When reset by 'R' input, it is said to be storing binary '0' ('Q' output = Low).


Like any other flip flop, the 'D' flip flop has two outputs that determine whether it is storing a binary '1' or a binary '0'. It also has two inputs. These are called 'D' and 'T' and work differently. The data or bit to be stored (which can be either a binary '0' or '1') is applied to the 'D' input. The 'T' input line controls the flip flop. It is used to determine whether the input data at 'D' is to be recognized or to be ignored. If the 'T' input is High, the data on the 'D' line get stored in the flip flop. If the 'T' line is low the 'D' input line data has no effect and the bit stored perviously is retained.

The 'JK' flip flop is the most versatile binary storage element. It can perform all the function of 'R', 'S' and 'D' flip flops plus it can do several other things. An integrated circuit 'JK' flip flop is really two 'RS' flip flops in one. These are called Master and Slave. Both flip flops are controlled by a common clock pulse to the 'T' input. When the 'T' line goes High, cutting off the Slave. At the same time data on 'J' and 'K' inputs is passed on to the Master for storage. When 'T' line goes Low, cutting off the Master from the input circuits. At the same time gates 'C' and 'D' are enabled and data stored in the Master is transferred to the Slave. This technique provides a complete isolation between the inputs and outputs. The integrated circuit IC 7476 contains two identical JK flip-flops which are completely independent except for a common power supply input connection.

PROCEDURE

Verification of 'RS' Flip Flop :

1. Connect the 4 logic inputs to 'Preset (PR)', Clear (CR)', 'S' & 'R' input of the Flip-Flop as shown in Fig. No. (1) through patchcords. Also connect 'Q' & 'Q' outputs to output indicators.
2. Connect 1Hz clock output to 'Clock (CK)' input of the flip flop.
3. Switch ON the instrument using ON/ OFF toggle switch provided on the front panel.
4. Vary the Truth Table No. (1) for


TRUTH TABLE No. (1) 'RS' FLIP FLOP

INPUTS					OUTPUTS	
PRESET (PR)	CLEAR (CR)	CLOCK (CK)	S	R	Q	$\bar{Q}$
L	H	X	X	X	H	L
H	L	X	X	X	L	H
L	L	X	X	X	H	H
H	H	P	L	L	Q	$\bar{Q}$
H	H	P	H	L	H	L
H	H	P	L	H	L	H
H	H	P	H	H	Not To Be Used	

various sets of input combinations.

Verification of 'D' Type Flip Flop :

1. Connect the output of NOT Gate to "R" input through patchcord as shown in Fig. No. (2). Connect 3 logic inputs to 'Preset (PR)', Clear (CR)' & 'D' input of the Flip-Flop as shown in Fig. No. (2) through patchcords. Also connect 'Q' & $\overline{Q}$ outputs to output indicators.
2. Connect 1Hz clock output to 'Clock (CK)' input of the flip flop.
3. Switch ON the instrument using ON/OFF toggle switch provided on the front panel.
4. Verify the Truth Table No. (2) for various sets of input combinations.


'D' TYPE FLIP FLOP
FIG. No. (2)

TRUTH TABLE No. (2) 'D' FLIP FLOP

INPUTS				OUTPUTS	
PRESET (PR)	CLEAR (CR)	CLOCK (CK)	D	Q	$\overline{Q}$
L	H	X	X	H	L
H	L	X	X	L	H
L	L	X	X	H	H
H	H	P	H	H	L
H	H	P	L	L	H

Verification of 'JK' Flip Flop :

1. Connect the 4 logic inputs to 'Preset (PR)', Clear (CR)', 'J' & 'K' input of the Flip-Flop as shown in Fig. No. (3) through patchcords. Also connect 'Q' & $\overline{Q}$ outputs to output indicators.

'JK' FLIP FLOP
FIG. No. (3)

TRUTH TABLE No. (3) 'JK' FLIP FLOP

INPUTS					OUTPUTS	
PRESET (PR)	CLEAR (CR)	CLOCK (CK)	J	K	Q	$\overline{Q}$
L	H	X	X	X	H	L
H	L	X	X	X	L	H
L	L	X	X	X	H	H
H	H	P	L	L	Q_0	$\overline{Q_0}$
H	H	P	H	L	H	L
H	H	P	L	H	L	H
H	H	P	H	H	TOGGLE	

2. Connect 1Hz clock output to 'Clock (CK)' input of the flip flop.

3. Switch ON the instrument using ON/ OFF toggle switch provided on the front panel.
4. Verify the Truth Table No. (3) for various sets of input combinations.

Verification of 'T' Type Flip Flop :

1. Short the "J" & "K" input of the IC 7476 to form "T" input. Also connect three logic inputs to 'Preset (PR)', Clear (CR)' & 'T' inputs of the Flip-Flop (To obtain "T" input short the 'J' & 'K' inputs). Also connect 'Q' & 'Q' outputs to output indicators.

2. Connect 1Hz clock output to 'Clock (CK)' input of flip flop.

3. Switch ON the instrument using ON/ OFF toggle switch provided on the front panel.

4. Verify the Truth Table No. (4) for various sets of input combinations.

'T' TYPE FLIP FLOP		
INPUT	OUTPUT	
	Q	Q
T	Q	Q
1	TOGGLE	

TRUTH TABLE No. (4)


Verification of 'D' Type Flip Flop :

1. Connect 3 logic inputs to 'Preset (PR)', Clear (CR)' & 'D' input of the Flip-Flop as shown in Fig. No. (4) through patchcords. Also connect 'Q' & 'Q' outputs to output indicators.

2. Connect 1Hz clock output to 'Clock (CK)' input of the flip flop.

3. Switch ON the instrument using ON/ OFF toggle switch provided on the front panel.

4. Verify the Truth Table No. (5) for various sets of input combinations.

'D' TYPE FLIP FLOP
FIG. No. (4)

TRUTH TABLE No. (5) 'D' FLIP FLOP

INPUTS				OUTPUTS	
PRESET (PR)	CLEAR (CR)	CLOCK (CK)	D	Q	Q
L	H	X	X	H	L
H	L	X	X	L	H
L	L	X	X	H	H
H	H	P	H	H	L
H	H	P	L	L	H

EXPERIMENT NO: 10

AIM: TO DESIGN & VERIFY OPERATION OF HALF ADDER & FULL ADDER.

APPARATUS REQUIRED: Power supply, IC's , Digital Trainer, Connecting leads .

BRIEF THEORY: We are familiar with ALU, which performs all arithmetic and logic operation but ALU doesn't perform/ process decimal no's. They process binary no's.

Half Adder: It is a logic circuit that adds two bits. It produces the O/P, sum & carry. The Boolean equation for sum & carry are

$$\text{SUM} = A + B$$

$$\text{CARRY} = A \cdot B$$

Therefore, sum produces 1 when A&B are different and carry is 1 when A&B are

1. Application of Half adder is limited.


Full Adder: It is a logic circuit that can add three bits. It produces two O/P sum & carry. The Boolean Equation for sum & carry are

$$\text{SUM} = A + B + C$$


$$\text{CARRY} = A \cdot B + (A + B) \cdot C$$

Therefore, sum produces one when I/P is containing odd no's of one & carry is one when there are two or more one in I/P.

CIRCUIT DAIGRAM

HALF ADDER

FULL ADDER

PROCEDURE:

- (a) Connect the ckt. as shown in fig. For half adder.
- (b) Apply diff. Combination of inputs to the I/P terminal.
- (c) Note O/P for Half adder.
- (d) Repeat procedure for Full wave.
- (e) The result should be in accordance with truth table.

OBSERVATION TABLE:

HALF ADDER:

INPUTS		OUTPUT	
A	B	S	C
0	0	0	0
0	1	1	0
1	0	1	0
1	1	0	1

FULL ADDER:

INPUTS			OUTPUTS	
A	B	C	S	CARRY
0	0	0	0	0
0	0	1	1	0
0	1	0	1	0
0	1	1	0	1
1	0	0	1	0
1	0	1	0	1
1	1	0	0	1
1	1	1	1	1

RESULT: The Half Adder & Full Adder ckts. are verified.

PRECAUTIONS:

- 1) Make the connections according to the IC pin diagram.
- 2) The connections should be tight.