Manual on Water and Wastewater Analysis

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INDEX

Experiment No.	Experiment Name	Date	Sign.
1	To determine pH of a given water/wastewater sample using pH electrode.		
2	To determine suspended, dissolved and total solids in a given water/wastewater sample by gravimetric method.		
3	To determine inorganic and organic solids in a given water/wastewater sample using gravimetric method.		
4	To determine acidity in a given water/wastewater sample.		
5	To determine alkalinity in a given water/wastewater sample.		
6	To calibrate nephelometer and determine turbidity in a given water sample.		
7	To determine chloride ion concentration in a given water /wastewater sample.		
8	To determine hardness in a given water sample by EDTA method.		
9	To determine optimum dose of alum for turbidity removal in a sample of water using jar test.		
10	To determine dissolved oxygen in a given water sample using modified winkler method.		
11	To determine biochemical oxygen demand of a given wastewater sample.		
12	To determine chemical oxygen demand of a given wastewater sample.		

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Aim: To calibrate pH meter and measurement of pH of given water sample.

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Learning Outcomes : The student will learn the operating Experimental Procedure required to calibrate a pH meter and will gain the experience for measuring pH of unknown water samples.

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Theory: Measurement of pH is one of the most important and frequently used tests in water chemistry. The pH of a solution is measured as negative logarithm of hydrogen ion concentration (Eq.1). $\text{pH} = -\log [\text{H}^+]$ Eq. (1) Where $[\text{H}^+]$ is the concentration (or activity) of hydrogen ion (or proton) in moles per liter (M). At a given temperature, the intensity of the acidic or basic character of a solution is indicated by pH or hydrogen ion concentration. pH values from 0 to 7 are diminishing acidic, 7 to 14 increasingly alkaline and 7 is neutral. The basic Brief Theory of electrometric pH measurement is determination of the activity of the hydrogen ions by potentiometric measurement using a standard hydrogen electrode and a reference electrode. Majority of modern pH electrodes are combination pH electrodes rather than a set of two separate electrodes. The combination pH electrode is virtually a tube within a tube arrangement. The inner tube housing the pH indicator electrode (pH sensing membrane, Ag/AgCl reference electrode, and HCl) and the outer one housing the reference electrode (Ag/AgCl) and its salt bridge. The pH sensing component of the indicator electrode is a glass bulb. When a pH probe is immersed in a sample solution, an electrical potential is developed inside and another

electrical potential is developed outside, the difference in the potential is measured and is given as the pH of the sample.

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Calibration of pH: Meter To be certain of accurate and reliable measurements, you need to perform pH meter calibration. This is generally done by measuring different buffer solutions with standardized, well-defined values, and then adjusting the pH meter based on any deviations from the buffer's known pH value. Buffers are solutions that resist changes in pH when small amounts of a strong acid or a base are added. Because they maintain a constant pH, buffers are good reference solutions. Calibration should be performed periodically to maintain the pH data quality.

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Apparatus: pH meter, Glassware (Beaker), tissue papers

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Reagents Buffer Solutions

1. pH 4: It is prepared by dissolving buffer tablet of pH 4 in a 100 mL volumetric flask and diluting it up to the mark with deionized water. Alternatively, it may be prepared by dissolving 1.021 g of potassium hydrogen phthalate ($\text{KHC}_8\text{H}_4\text{O}_4$) in fresh 100 mL deionized water.
2. pH 7: Dissolve a buffer tablet of pH 7 in deionized water in a 100 mL volumetric flask and make up to the mark. Alternatively, dissolve 0.340 g of potassium hydrogen phosphate (KH_2PO_4) and 0.355 g disodium hydrogen phosphate

(Na₂HPO₄) in distilled water in a 100 mL volumetric flask and make up the solution to the mark.

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Procedure

Calibration of pH meter

1. Plug in the meter and let it warm up for five to ten minutes.
2. Put the selector knob to pH mode.
3. Rinse the electrode with de-ionized water and blot dry using a piece of tissue.
4. Allow all the buffers to reach the same temperature, since pH readings are temperature dependent. If the buffers are not at 25 °C, temperature compensation is recommended. Measure the temperature of the buffers using thermometer and manually adjust the “temperature compensation knob” to this temperature. Where, ATC probe is available the system automatically transmits the temperature of the buffers to the meter.
5. Place the electrode in the solution of pH 7 buffer. It should be ascertained that the glass electrode membrane is completely immersed in the solution. Electrode should not touch the sides or the bottom of the beaker. Allow the pH value to stabilize and adjust 'set buffer knob' in a manner that the reading matches with buffer pH. Remove the electrode from the buffer.
6. Rinse the electrode with de-ionized water and blot dry using a piece of tissue.
7. Place the electrode in the solution of next buffer i.e. pH 4. Follow the step no. 4 until pH value reaches pH 4. Remove the electrode from the buffer. Similarly, other pH standard such as 10.1 can be used, depending on the expected sample pH range. Generally, for groundwater sample pH 7.0 and 10.0 should be used for calibration purpose.
8. Rinse the electrode with de-ionized water and blot dry using a piece of tissue.
9. Now, the instrument is ready to measure pH of unknown sample.

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pH Measurement of samples

1. Make sure that the meter is set to the pH mode.
2. Place the electrode in the sample (30-50 mL) to be tested.
3. The pH of the solution appears in the display. Allow the display to stabilize before taking your reading.
4. Note down the pH and temperature reading.
5. Rinse the pH electrode and place it back in the storage solution.

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Precautions

1. Make sure to store a pH electrode in its wetting cap containing electrode fill solution (3 M KCl, purchased or prepared by dissolving 22.37 g KCl into 100 mL DI water). Do not store electrode in DI water.
2. Never touch the membrane of the glass electrode with anything else except soft tissue paper since it is fragile and is easily ruined if scratched or bumped.
3. Always bring the standards and the sample to the same temperature before measurement (preferably room temperature).

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Questions:

Q1. Write down the environmental significance of pH measurement of water samples.

Q2. One solution has a pH of 4.0 and another has pH of 6.0. What is (a) the hydrogen ion activity and (b) the hydroxide ion activity in each of the solutions?

Q3. A 50 % decrease in hydrogen ion activity represents how much of an increase in pH units?

Experiment No. – 2 & 3

2 & 3 Total Solids, Total Suspended Solids, Total Dissolved Solids, organic & inorganic solids

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Aim: To determine Total Solids, Total Suspended Solids, Total Dissolved Solids, organic & inorganic solids in water or wastewater sample through Gravimetric Method

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Background:

Solids analysis provides one of the fundamental measurements used for control of the activated sludge process and for the regulation of wastewater discharges.

Gravimetric analysis is based on the determination of constituents or categories of materials by measurement of their weight. The experiment illustrates the principles of weighing and demonstrates separation and categorization techniques used to define the various types of solids in waters and wastewaters. These techniques involve three analytic operations in addition to weighing. These are: filtration, evaporation, and combustion. *Filtration* is used to separate suspended or particulate (non-filterable) fraction from dissolved or soluble (filterable) fractions. *Evaporation* separates water from material dissolved or suspended in it. *Combustion* differentiates between organic and inorganic matter. Organic matter will be destroyed completely by burning at 550°C for 30 min.

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Procedure:

Samples:

1. Weigh filters (mass=B g).
2. Filter samples (50 ml).
3. Run each sample in duplicates. You will have a total of 4 samples.
4. Oven dry at 103°C for 30 min (please note that standard methods recommend 1 hour). At this stage all water will be evaporated and only suspended solids will be retained on filter. Weigh filters now (mass =A g).
5. **Calculate concentration of total suspended solids**
$$\text{mg total suspended solids} = 1000 * (A - B) / (\text{sample volume in mL}) \quad (1)$$
6. Weigh crucibles (mass=C) and then weight crucible and filter (total mass=E).
7. Put crucibles with filter paper in muffle furnace and Ignite at 550°C for 15 min. At this stage all volatile components of solids will be volatilized and only fixed inert solid materials will be left in crucible. Weigh crucible after proper cooling (mass=D g).
8. **Calculate concentration of volatile suspended solids (F) :**
$$\text{mg volatile suspended solids/L} = 1000 * (E - D) / (\text{sample volume in mL}) \quad (2)$$
9. **Calculate concentration of fixed suspended solids (G) :**
$$\text{mg volatile suspended solids/L} = 1000 * (D - C) / (\text{sample volume in mL})$$

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Acidity (अम्लता)

Objective :

To determine the acidity and alkalinity of a given water sample.

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Outcome : Students will learn the titration method to measure the acidity of water and wastewater samples.

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Theory: All waters having pH lower than approximately 8.3 should contain acidity. Acidity is the capacity of water to neutralize a base (OH). Sources of acidity in natural water include CO₂ (air, bacterial degradation), H₂PO₄, H₂S, protein, fatty acids, and salt of trivalent metals such as hydrated Fe³⁺ and Al³⁺. Acids contribute to corrosiveness and influence chemical reaction rates, chemical speciation and biological processes. Acidity of water is its quantitative capacity to react with a strong base to a designated pH. In polluted water, acidity may be caused by free mineral acid (H₂SO₄, HCl) from metallurgical industry, acid rain, and organic acid mine drainage, waste. To measure acidity, an alkaline solution (NaOH) and two indicators are used. The measured value may vary significantly with the end point pH used in the determination. The first indicator (methyl orange) changes colour from red to yellow at approximately pH 3.7. This corresponds to mineral acidity or methyl orange acidity. If NaOH is continued to titrate, the second indicator phenolphthalein changes colour from colourless to pink at pH 8.3. Titration to the phenolphthalein end point of pH 8.3 measures acidity from both mineral acids and weak acids. This total acidity is also termed phenolphthalein acidity. To calculate mineral acidity or total acidity from the titration data, the amount of NaOH must be converted to the equivalent amount of CaCO₃ (in mg/L) in water. Since the molecular weight (MW) of CaCO₃ is 100 g/mol or 1*10⁵ mg/mol, its equivalent weight is ½ MW or 5*10⁴. This is used as the conversion factor so that the acidity can be reported as mg CaCO₃/L. Note: As far as water analysis is concerned, acidity test does not bear significant importance because methyl orange acidity invariably remains absent in the raw water and even phenolphthalein acidity (that too principally due to the excessive-prevalence of dissolved carbon dioxide and carbonic acids) normally does not exist to a significant extent in the raw water.

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Apparatus Titration Stand with Burette, Conical Flask, Beaker, Measuring Cylinder, Pipette, Funnel, Wash Bottle etc.

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Reagents

- 0.02 N NaOH Solution: Dissolve 0.8 grams of NaOH in 1 L of deionized water.
- Phenolphthalein Indicator: Dissolve 0.5 gm of Phenolphthalein Indicator in 500ml 95% of ethyl alcohol. Add 500 mL distilled water.
- Methyl orange Indicator: Dissolve 0.5 g of Methyl orange Indicator & dilute to 1000 ml with distilled water.

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Experimental Procedure

Methyl orange acidity

1. Take suitable volume of sample (25, 50 or 100 mL) in a conical flask (V)
2. Add 2 drops of methyl orange and see the colour. If colour turns yellow, methyl orange acidity is absent. If colour turns orange, titrate with 0.02 N NaOH till colour changes to faint orange. Characteristic of pH 4.3-4.4. Note the volume of NaOH required (A).

Phenolphthalein Acidity

1. Take suitable volume of sample (25, 50 or 100 mL) in a conical flask (V).
2. Add 2-3 drops of phenolphthalein indicator. If the sample turns pink colour it means that, phenolphthalein acidity is absent and stop the experiment.
3. If there is no colour change, titrate with 0.02N Standard NaOH solution till faint pink colour appears indicating pH 8.3. The end point is colourless to faint pink colour. Note additional volume of NaOH required (B).

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Indicators: Phenolphthalein and Methyl Orange Indicator.

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End Point: a) Yellow to Faint Orange b) Colourless to Faint Pink Colour

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Observation Table:

S.No.	Volume of Sample (ml)	Initial Burette Reading (ml)	Final Burette Reading (ml)	Volume of NaOH (ml) used

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Mineral Acidity as CaCO₃ in mg/L =

[Volume of NaOH (A) × N × 50 × 1000] / Volume of sample]

Phenolphthalein Acidity as CaCO₃ in mg/L =

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Experiment No. - 5

Alkalinity (क्षारीयता)

Objective: Calculate hydroxide, carbonate and bicarbonate alkalinity.

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Theory:

In natural water, alkalinity is due to H_2CO_3 (CO_2 dissolved in water). Alkalinity from other materials is insignificant and may be ignored. In polluted waters, minor amounts of NH_3 and salts of weak acids such as borate, silicate, and phosphate (i.e., the conjugate bases of HBO_3 , H_4SiO_4 , H_3PO_4), and organic acids will contribute to alkalinity. Alkalinity of water is its quantitative capacity to neutralize a strong acid to a designated pH. In other words, its ability to maintain a relatively constant pH. The possibility to maintain constant pH is due to the hydroxyl, carbonate and bicarbonate ions present in water. The ability of natural water to act as a buffer is controlled in part by the amount of calcium and carbonate ions in solution. Carbonate and calcium ion come from calcium carbonate or limestone. So, water that comes in contact with limestone will contain high levels of both Ca^{++} and CO_3^{2-} ions and have elevated hardness and alkalinity. Alkalinity is significant in many uses and in treatment of natural acid and wastewater. It is significant in determining suitability of water for drinking as well as irrigation purpose. The titration Experimental Procedure can distinguish two types of alkalinity, that is, phenolphthalein alkalinity and total alkalinity (methyl orange alkalinity). The phenolphthalein alkalinity is the acid-neutralizing power of hydroxide and carbonate ions (OH^- and CO_3^{2-}) present in the water sample, while total alkalinity represents all the bases in it (OH^- , CO_3^{2-} and HCO_3^-). The titration Experimental Procedure for alkalinity is exactly the opposite of the Experimental Procedure for acidity. Two same indicators are used but in the reverse order. The Experimental Procedure utilizes the first indicator, phenolphthalein, to signal the end points (pH 8.3) in titrating OH^- and CO_3^{2-} . Further titration with H_2SO_4 reaches the end point of the second indicator (methyl orange) at around pH 4.5, which measures the total alkalinity. Note that the pH for the colour change is not at a fixed value, but rather a range, that is, pH 3.1–4.4 and pH 8.0–9.6 for methyl orange and phenolphthalein, respectively.

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Reagents:

- 0.02N H_2SO_4 : 3 mL. Conc. H_2SO_4 in 1 L deionized water to get 0.01 N H_2SO_4 . Phenolphthalein Indicator: Dissolve 0.5 gm of Phenolphthalein Indicator in 500 ml 95% of ethyl alcohol. Add 500 mL distilled water.
- Methyl orange Indicator: Dissolve 0.5 g of Methyl orange Indicator & dilute to 1000 ml with distilled water.

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Experimental Procedure:

1. Take suitable volume of sample (25, 50 or 100 mL) in a conical flask.
2. Add 2-3 drops of phenolphthalein indicator.

3. If the sample turns pink in colour, then titrate with 0.02 N H₂SO₄, till the pink colour disappears. Note the volume of acid required as (A).
4. If the pink colour does not exist it indicates that, phenolphthalein Alkalinity is absent. Then continue the titration with Methyl Orange indicator. Methyl Orange Alkalinity
5. Add 2-3 drops of methyl orange to the same flask; continue titration till colour becomes orange. Note down volume of H₂SO₄ consumed as 'B'.

S.No.	Volume of Sample (ml)	Initial Burette Reading (ml)	Final Burette Reading (ml)	Volume of H ₂ SO ₄ (ml) used

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Calculations:

Phenolphthalein Alkalinity as CaCO₃ in mg/L =

$$[(\text{Volume of NaOH(A)} \times N \times 50 \times 1000) / \text{Volume of sample}]$$

Methyl Orange Alkalinity as CaCO₃ in mg/L =

$$[(\text{Volume of NaOH (B)} \times N \times 50 \times 1000) / \text{Volume of sample}]$$

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Turbidity (टर्बिडिटी)

Aim

To calibrate and standardize Nephelometer and to determine the turbidity of a given water sample.

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Outcome

Students will learn the use and working Brief Theory of Nephelometer and will learn how to calculate turbidity of water samples.

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Theory

The Nephelometry method is based upon a comparison of the intensity of light scattered by the sample under defined conditions with the intensity of light scattered by a standard reference suspension. The higher the intensity of scattered light, the higher is the turbidity. Readings of the Nephelometer are given in Nephelometric Turbidity Unit (NTU). A primary standard suspension is used to calibrate the instrument. A secondary standard suspension is used as a daily calibration check and is monitored periodically for deterioration using one of the primary standards.

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Reagents

- Dissolve 1.000 g hydrazine sulfate $[(\text{NH}_2)_2 \cdot \text{H}_2\text{SO}_4]$ in filtered water and dilute to 100 mL in a volumetric flask.

- Dissolve 10.00 g hexamethylenetetramine $[(CH_2)_6N_4]$ in filtered water and dilute to 100 mL in a volumetric flask.
- Mix 5.0 mL of hydrazine sulfate and 5.0 mL of hexamethylenetetramine solutions in a 100-mL volumetric flask and let stand 24 hours at $25 \pm 3^\circ C$; dilute to the mark and mix. To prepare 500 mL of 400 NTU standard, mix 25 mL of the reagent solutions in a 500-mL flask, dilute to the mark, and mix.
- For a 40 NTU standard, dilute 10.00 mL of the 400 NTU stock suspension to 100 mL with turbidity-free water (sample or deionized water passed through a filter media of 0.2 mm).

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Experimental Procedure:

(A) Calibration of Nephelometer

The calibration of Nephelometer can be done using the following steps:

1. Switch the turbidimeter on and allow it to warm up.
2. Check instrument focus: insert template in the cell holder. The lamp image should just fill the inside circle.
3. Draw out the sample cell from the Nephelometer and hold it by the rim not beneath the lip.
4. Pour the deionized water into the sample cell and wipe out the exterior with a soft dry tissue paper.
5. Using Set Zero dial calibrate the instrument to Zero.
6. Next, pour the standard solution into the sample cell and wipe out the exterior with a soft dry tissue paper. Place the cell in its holder with correct orientation.

7. Select the desired NTU range dial of the Nephelometer and rotate the calibrate dial of the Nephelometer to match the digital reading exactly with the preset range.
8. The instrument is now calibrated.
9. The Sample cell can now be filled with unknown sample to measure the NTU. Calculation

The NTU of Diluted sample can be calculated using the following formula:

$$NTU = A \times (B+C) / C$$

Where,

A = NTU found in diluted sample,

B = volume of dilution water, in milliliters, and

C = sample volume taken for dilution, in milliliters.

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Results:

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Home Work

- Q1. Write down the environmental significance of measuring turbidity of water samples.
Q2. Why turbidity is important in "filtration" and "disinfection" processes?
Q3. Discuss why turbidity in general cannot be correlated with weight concentration of suspended matter in water samples?

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Experiment No. - 7

Chloride(क्लोराइड)

Aim: To determine chloride ion concentration in a water sample.

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

Outcomes: Students will learn to measure the chloride concentration of water samples.

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Theory :


Chlorides are present in water usually as NaCl, MgCl₂ and CaCl₂. Although chlorides are not harmful as such, their concentrations over 250 ppm impart a peculiar taste to the water thus rendering the water unacceptable for drinking purposes. By Mohr's method, chloride ions in a water sample (neutral or slightly alkaline) can be determined by titrating it against standard silver nitrate (AgNO₃) solution using potassium chromate (K₂CrO₄) as an indicator. The pH should be in between 7-8. At higher pH, silver ions are precipitated as silver hydroxide. At lower pH, potassium chromate indicator is converted to potassium dichromate (K₂Cr₂O₇). Mohr's method is based on the precipitation titration in which silver nitrate solution is released from the burette to the water sample which contains chloride ions and indicator. The silver ions (from silver nitrate solution) react with chloride ions (from water sample) and chromate ions (from indicator) to form white precipitate of silver chloride and red precipitate of silver chromate. The following reactions depict the detection of chlorine.

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Apparatus

Burette, conical flask, pipette, measuring cylinder.

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Reagents

- Standard silver nitrate solution (N/50),
- Indicator potassium chromate solution.

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Calibration:-

Transfer 25 ml of the distilled water in a conical flask and add 3-4 drops of indicator potassium chromate solution. Slowly add standard silver nitrate solution from the burette and shake the solution well. At the end point, light yellow colour starts changing to red colour. The titration is repeated until a concordant volume is obtained. The blank correction for the indicator should be subtracted from the volume of the titrant obtained after titrating the with the unknown water sample solution. Let us consider the titrated volume as (V1).

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Procedure

- Take 50 ml sample in a conical flask. Measure sample pH (to check for acidity or basicity of the sample, this step is optional).
- Add 1.0 ml of freshly prepared indicator solution of $K_2Cr_2O_4$.
- Titrate with standard silver nitrate solution until the brick red precipitate appears. This stage marks the end point of titration and note down volume of titrant used. Note the difference of volume as (V2). Find the actual volume (V) of titrate used using: $V = V_2 - V_1$ equation.

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Observation

S.No.	Volume of Sample (ml)	Initial Burette Reading (ml)	Final Burette Reading (ml)	Volume of $AgNO_3$ (ml) used

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Where, A = volume of titrant used, N is normality of silver nitrate (here we used N/50 or 0.02 N), and V_{sample} is volume of sample used (mL).

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$$\text{क्लोराइड आयन सांद्रता (मिग्रा/लीटर)} = (A \times N \times 35.45) \times 1000 / V_{\text{sample}}$$

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

Aim: To determine the hardness of a water sample

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

Learning Outcomes Students will learn to measure the hardness of the water samples. Brief Theory All-natural waters have salts dissolved in them. These problems are caused mainly by the presence of calcium and magnesium ions in the water. Water hardness is caused by multivalent metallic cations, primarily the divalent Ca^{2+} and Mg^{2+} . Anions and monovalent cations such as Na^+ and K^+ do not contribute to hardness. Hard water adversely affects the suitability of water for uses in domestic and industrial purposes. You may have experienced this problem if you have taken a bath with water that contained too many calcium and magnesium salts. Eq. 1 and 2 present the respective hardness.

Although hardness is caused by cation, it may also be discussed in terms of carbonate (temporary) and non-carbonate (permanent) hardness. Carbonate hardness refers to the amount of carbonates and bicarbonates in solution that can be removed or precipitated by boiling. This type of hardness is responsible for the deposition of scale in hot water pipes and kettles. When total hardness is numerically greater than that of total alkalinity expressed as CaCO_3 , the amount of hardness equivalent to alkalinity is called carbonate hardness¹. When the hardness is numerically equal to less than total alkalinity, all hardness is carbonate hardness. The amount of hardness in excess of total alkalinity expressed as CaCO_3 is non- carbonate hardness. Non carbonate hardness is caused by the association of the hardness- causing cation with sulphate, chloride or nitrate and is referred to as “permanent hardness”. This type of hardness cannot be removed by boiling. Hardness is determined by the EDTA method in alkaline condition; EDTA and its sodium salts form a soluble chelated complex with certain metal ions as shown in following reaction;

The reaction between calcium ions and EDTA only occurs at a high pH; therefore, a solution that has a constant pH of 10 (called a buffer) will be added to the water sample. Eriochrome black T, a metal ion indicator, will be used to visualize the endpoint. Eriochrome Black T (EBT) (blue colour solution) serves as an excellent indicator to show when all hardness ions have been consumed. When small amount of EBT is added to hard water with pH>10, it combines with Ca²⁺ and Mg²⁺ ions to form weak complex ions (wine-red colour solution). A small amount of indicator is added to the solution containing the Ca²⁺ forming a wine-red complex.

The end point will be indicated when the original red solution turns to blue indicating that the EDTA has reacted with all the calcium ions in the water sample. In this experiment you will measure calcium carbonate (CaCO₃) in ppm. You will be able determine the hardness of your water sample by measuring the volume of your water sample and the volume (in mL) of the EDTA solution used to react with all of the calcium in the water sample.

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Experimental Procedure

1. Take 50 mL sample in a conical flask
2. Using a dropper add few drops of the pH 10 buffer solution to the conical flask.
3. Using spatula, add a pinch of indicator Eriochrome black T (in the form of a powder) to the flask. The colour of the solution will be red.
4. Clean and fill burette with the EDTA solution and record the initial volume.
5. Titrate the solution with EDTA till the colour of solution will changed from red to blue, record the final volume of burette and put in the given formula for calculating the hardness.

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Observation Table:

S.No.	Volume of Sample (ml)	Initial Burette Reading (ml)	Final Burette Reading (ml)	Volume of NaOH used

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Calculations

Hardness as per CaCO₃ (mg/L) =

{mL of EDTA used in titration × 1000 mg/L} / mL of Water sample

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Outcome Students will learn the use and working Brief Theory of Nephelometer and will learn how to calculate turbidity of water samples.

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Theory The Nephelometry method is based upon a comparison of the intensity of light scattered by the sample under defined conditions with the intensity of light scattered by a standard reference suspension. The higher the intensity of scattered light, the higher is the turbidity. Readings of the Nephelometer are given in Nephelometric Turbidity Unit (NTU). A primary standard suspension is used to calibrate the instrument. A secondary standard suspension is used as a daily calibration check and is monitored periodically for deterioration using one of the primary standards.

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Reagents

- Dissolve 1.000 g hydrazine sulfate $[(\text{NH}_2)_2 \cdot \text{H}_2\text{SO}_4]$ in filtered water and dilute to 100 mL in a volumetric flask.
- Dissolve 10.00 g hexamethylenetetramine $[(\text{CH}_2)_6\text{N}_4]$ in filtered water and dilute to 100 mL in a volumetric flask.
- Mix 5.0 mL of hydrazine sulfate and 5.0 mL of hexamethylenetetramine solutions in a 100-mL volumetric flask and let stand 24 hours at $25 \pm 3^\circ\text{C}$; dilute to the mark and mix. To prepare 500 mL of 400 NTU standard, mix 25 mL of the reagent solutions in a 500 mL flask, dilute to the mark, and mix.
- For a 40 NTU standard, dilute 10.00 mL of the 400 NTU stock suspension to 100 mL with turbidity-free water (sample or deionized water passed through a filter media of ≤ 0.2 mm).

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- Switch the turbidimeter on and allow it to warm up.
- Check instrument focus: insert template in the cell holder. The lamp image should just fill the inside circle.
- Draw out the sample cell from the Nephelometer and hold it by the rim not beneath the lip.
- Pour the deionized water into the sample cell and wipe out the exterior with a soft dry tissue paper.
- Using Set Zero dial calibrate the instrument to Zero.
- Next, pour the standard solution into the sample cell and wipe out the exterior with a soft dry tissue paper. Place the cell in its holder with correct orientation.
- Select the desired NTU range dial of the Nephelometer and rotate the calibrate dial of the Nephelometer to match the digital reading exactly with the preset range.
- The instrument is now calibrated.
- The Sample cell can now be filled with unknown sample to measure the NTU.

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Experiment No. – 10

DO

Aim/Objective - To estimate Dissolved Oxygen (DO) level in the given water sample by Winkler's Method.

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Theory –

- a) Dissolved Oxygen or DO, as the name suggests is the amount of free oxygen dissolved in the water which can be used by the aquatic life.
- b) The solubility of oxygen in water depends on temperature and other factors. The presence of oxygen is essential for the survival of aquatic life in water. A rapid fall of DO in river water is one of the first indications of pollution.
- c) Two methods are commonly used to determine to DO concentration:
 - i. The iodometric method which is a titration-based method and depends on the oxidizing property of DO
 - ii. The membrane electrode procedure, which works based on the rate of diffusion of molecular oxygen across a membrane.

Winkler's Modified Method:

- a. The Winkler's Method is the technique used to measure DO in freshwater structures. It is used as an indicator of the health of a water body, where higher DO concentrations are associated with high production and little contamination.

- b. This test is performed on-site, as delays between sample collections and testing may result in a variation in oxygen content.
- c. In this method, the oxygen present in the water sample oxidizes the divalent manganous ion to its higher valency, which precipitates as a brown hydrated oxide after addition of sodium hydroxide (NaOH) and potassium iodide (KI).
- d. Upon acidification, manganese reverts to divalent state and liberates iodine from KI equivalent to DO content in the sample. The liberated iodine is titrated against sodium thiosulphate (0.025N), using freshly prepared 2% starch solution as indicator. If the oxygen is absent in the sample, the MnSO_4 react with the alkali to form white precipitate Mn(OH)_2 .
- e. When MnSO_4 and alkali-iodide reagent ($\text{NaOH} + \text{KI}$) are added to a sample in the absence of oxygen, a pure white precipitate is formed. $\text{Mn}^{2+} + 2\text{OH}^- \rightarrow \text{Mn(OH)}_2$ (white precipitate). Mn^{2+} is oxidized to Mn^{4+} and precipitates brown hydrated oxide in the presence of oxygen. $\text{Mn}^{2+} + 2\text{OH}^- + 0.5\text{O}_2 \rightarrow \text{MnO}_2$ (brown hydrated precipitate) + H_2O . The oxidation of Mn^{2+} to MnO_2 is called fixation of the oxygen, and occurs slowly at low temperature. $\text{Mn(OH)}_2 + 0.5\text{O}_2 \rightarrow \text{MnO}_2 + \text{H}_2\text{O}$
- f. After the floc has settled, sulfuric acid is added. MnO_2 is oxidized to produce I_2 under low pH conditions. I_2 is insoluble in water and forms complexes when there is an excess of iodide ions in solution, preventing iodine ions from escaping. $\text{MnO}_2 + 2\text{I}^- + 4\text{H}^+ \rightarrow \text{Mn}^{2+} + \text{I}_2 + 2\text{H}_2\text{O}$ and alongside $\text{I}_2 + \text{I}^- \leftrightarrow \text{I}_3^-$. The sample is now ready to be titrated with thiosulfate solution.

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Procedure –

a. Standardisation of Sodium Thiosulphate Solution –

- Fill the burette with the Sodium Thiosulphate Solution ($Na_2S_2O_3$) up to the 0 mark with the help of a funnel.
- Now in a conical flask we add around 10 mL of 0.025 N Potassium Dichromate ($K_2Cr_2O_7$) with the help of a graduated pipette, giving the solution a yellow colour.
- To this beaker and solution, we add around 2 mL of Hydrochloric Acid (HCl) and 1 spatula of Potassium Iodide (KI). This turns the colour of the solution to a brown-yellow colour.
- Now, we titrate this solution against the Sodium Thiosulphate Solution ($Na_2S_2O_3$) till a pale-yellow colour is achieved. Post this colour, we add around 2 mL of Starch to the beaker turning the beaker colour to blue. Finally, it is again titrated till the final colour becomes white.
- Finally, using the formula $N_1V_1 = N_2V_2$, we achieve the Normality of the Sodium Thiosulphate Solution ($Na_2S_2O_3$).

b. Determination of Dissolved Oxygen of Sample –

- Fill the BOD bottle with the test sample up to the brim of the bottle.
- Immediately after this we add around 2 mL of Manganese Sulphate ($MnSO_4$) solution to this bottle, just below the surface of the sample.

- iii. Beyond this we also add 2 mL of Alkali-Iodize-Azide reagent to this bottle similar to above step. This turns the solution to a golden-brown colour.
- iv. Mix this solution properly by inverting it several times and allow the precipitate to settle to the bottom by leaving it stagnant for a few minutes.
- v. Now carefully add 2mL of Concentrated Sulphuric Acid (H₂SO₄). Again, mix the solution by inverting the bottle several times and then finally transfer 200 mL of this solution to a beaker.
- vi. Titrate this solution against the Sodium Thiosulphate Solution (Na₂S₂O₃) till a pale-yellow colour is achieved. Post this colour, we add around 2 mL of Starch to the beaker turning the beaker colour to blue. Finally, it is again titrated till the final colour becomes white.
- vii. Finally, using the formula given below we can find the Dissolved Oxygen content in the solution.

$\text{Dissolved Oxygen (D.O.) as mg/l} = \frac{N_{\text{Na}_2\text{S}_2\text{O}_3} \times V_{\text{Na}_2\text{S}_2\text{O}_3} \times 8000}{\text{Volume of the sample}}$ $N_{\text{Na}_2\text{S}_2\text{O}_3}$ = Normal of Na₂S₂O₃ $V_{\text{Na}_2\text{S}_2\text{O}_3}$ = Volume of Na₂S₂O₃

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

b) Standardisation of Silver Nitrate –

- i. Initial burette reading =
- ii. Final burette reading =
- iii. Volume of Sodium Thiosulphate Solution ($\text{Na}_2\text{S}_2\text{O}_3$) used =
- iv. Volume of Potassium Dichromate =

Results –

- a. Normality of Sodium Thiosulphate Solution ($\text{Na}_2\text{S}_2\text{O}_3$) =
- b. Dissolved Oxygen Content =
- c. The acceptable limit of DO by BIS is 4 – 7 mg/L. DO of this sample is mg/L thus it is within safe limits and is safe for drinking purpose.

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Objective: To determine BOD value for determining biodegradability of solution.

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Background: The most widely used test indicating organic pollution of both wastewater and surface water is the 5-day BOD (BOD₅). This determination involves the measurement of the dissolved oxygen used by microorganisms in the biochemical oxidation of organic matter. BOD₅ is the total amount of oxygen consumed by microorganisms during the first five days of biodegradation. Oxygen demand is associated with the biodegradation of the carbonaceous portion of wastes and oxidation of nitrogen compounds such as ammonia.



The following equations simplify the process of biodegradation:

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Procedure:

Apparatus: Incubation bottle 300mL volume; Air compressor, 20°C incubator

Reagents for DO measurement:

- Manganese sulfate solution: Dissolve 480 g MnSO₄·4H₂O, 400 g MnSO₄·2H₂O or 364 g MnSO₄·H₂O in distilled water, filter, and dilute to 1L. The MnSO₄ solution should not give a color with starch when added to an acidified potassium iodide (KI) solution.

- b) Alkali-iodide-azide reagent
- c) Sulfuric acid: One mL is equivalent to ~ 3mL alkali-iodide-azide reagent.
- d) Starch solution: Dissolve 2 g laboratory-grade soluble starch and 0.2 g salicylic acid as preservative in 100 mL hot distilled water.
- e) Standard sodium thiosulfate titrant: Dissolve 6.205 g $\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$ in distilled water and add 1.5 mL 6N NaOH or 0.4 g solid NaOH and dilute to 1000 mL. Standardize with bi-iodate solution.
- f) Standard potassium bi-iodate solution (0.0021M): Dissolve 812.4 mg $\text{KH}(\text{IO}_3)$ in distilled water and dilute to 1000 mL.
- g) Standardization: Dissolve ~ 2 g KI, free from iodate in an Erlenmeyer flask with 100 to 150 mL distilled water; add 1 mL 6N H_2SO_4 or a few drops of conc. H_2SO_4 and 20.00 mL standard bi-iodate solution. Dilute to 200 mL and titrate liberated iodine with thiosulfate titrant, adding starch toward end of titration, when a pale straw color is reached. When the solution is of equal, 20.00 mL 0.025M $\text{Na}_2\text{S}_2\text{O}_3$ should be required. If not, adjust the $\text{Na}_2\text{S}_2\text{O}_3$ solution to 0.025M.

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Steps:

DO measurement:

1. Make dilution water by adding 2mL/L of following reagents in distilled water:
 - a. Phosphate buffer solution
 - b. Magnesium sulfate solution
 - c. Calcium chloride solution
 - d. Ferric chloride solution
 - e. Sodium Sulfite solution
2. For a given sample bottle, add 1 mL of alkali azide and then 1 mL manganous sulfate solution. Shake well the bottle and keep it open for 5 minutes to settle the precipitate. Add 2 mL concentrated H₂SO₄ and place the cap on the bottle. Shake well the bottle till all the precipitate is dissolved.
3. Take 203 mL of sample in conical flask and titrate with standard sodium thiosulfate solution (0.025N) till the colour changes from dark yellow to light yellow. Then add few drops of starch indicator and continue to titrate till the color of the solution becomes either colorless or changes to its original sample colour. Note down volume of 0.025N sodium thiosulfate consumed.
4. Calculate DO value of the sample. Remember that in 200 mL sample, 1 mL of sodium thiosulfate of 0.025N equals to 1 mg/L dissolved oxygen:
=>Dissolved oxygen (DO) (in mg/L) = mL of sodium thiosulfate (0.025N) consumed.

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BOD:

1. Prepare BOD dilutions. Use dilution water (it contains nutrients, the exact contents are described in Standard Methods): Blank (only dilution water); 5 mL sample in 300 mL BOD bottle, fill up with dilution water; 15 mL sample in 300 mL BOD bottle, fill up with dilution water; 20 mL sample in 300 mL BOD bottle, fill up with dilution water

2. Take 300 mL sample in BOD bottle. Prepare two sets of this sample. Keep one set for DO analysis for day 0 (i.e., Sample0Day) and another sample in BOD incubator for 5 days at 20° C (Sample5Day).

3. Measure DO in different samples at $t=0$.

4. Incubate samples in 20°C for 5 days.

5. Come back in the lab after 5 days and record dissolved oxygen.

6. Record data in following manner. From table

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Calculate 5-day BOD value of the sample at 20°C:

$$\text{t-day BOD} = [\text{DO}_t - \text{DO}_0] / (P) \quad (1)$$

where P= Dilution factor = 300mL/(sample volume in mL)

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Experiment No. -12 Chemical Oxygen Demand (COD)

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Objective: To determine COD value for determining organic strength of solution (Closed Reflux Method).

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Background:

Chemical oxygen demand (COD) is termed as the amount of a specific oxidizing agent that reacts with sample under controlled conditions and it is expressed as oxygen equivalence. This parameter indicates the extent of organic matter contamination of water and is always higher than the biochemical oxygen demand (BOD). It is used to indicate organic matter contamination and it helps in knowing overall organic load to the receiving body.

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Selection of Method

There are two methods for COD determination.

The first method: open reflux method is suitable for a wide range of wastes where large volume of sample is required (for samples with COD= 50 mg O₂/L).

In the second method: closed reflux methods, small quantities of metallic salt reagents are required and small quantities of hazardous waste is produced (for samples with COD= 5 to 50 mg O₂/L). In the closed reflux method, ampules and culture tubes with premeasured reagents are used and then samples is placed in the tube and COD is determined.

In this experiment, closed reflux method is used and samples with COD < 50 mg O₂/L are tested.

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

Reaction with dichromate solution of sample:

Potassium dichromate is a strong oxidizing agent and it can be used to prepare solution of exact normality.

$$\text{Here } d = (2n/3) + (a/6) - (b/3) - (c/2)$$

During experiment, excess dichromate concentration is determined by titrating it with ferrous ammonium sulfate (FAS). The reaction is given by:


$$\text{Here } d = (2n/3) + (a/6) - (b/3) - (c/2)$$

Ferriin (ferrous 1,10-phenanthroline sulfate):

It is used to indicate change in oxidation-reduction potential of the solution and it indicates the condition when all dichromate has been reduced by ferrous ion. It gives a very sharp brown color change which can be seen in spite of blue color generated by the Cr³⁺ ions formed on reduction of the dichromate.

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
Procedure:

Apparatus: Digestion vessels; block heater; microburet; ampule sealer. Borosilicate culture tubes (16mm*100 mm or 20 mm*150mm) with TFE lined-screw caps are used. The block heater is required to operate at $150 \pm 2^\circ\text{C}$ with holes to accommodated digestion vessels. **Do not use an oven because of the possibility of leaking samples generating corrosive and explosive atmosphere.**

Reagents:

a. Standard potassium dichromate digestion solution, 0.01667M: Add to about 500 mL distilled water 4.903 g $\text{K}_2\text{Cr}_2\text{O}_7$, primary standard grade, previously dried at 150°C for 2 h, 167 mL conc. H_2SO_4 , and 33.3 g HgSO_4 . Dissolve, cool to room temperature, and dilute to 1000 mL.

b. Sulfuric acid reagent:

c. Ferroun indicator solution: Dilute it by a factor of 5 as required. This indicator is used to indicate change in oxidation-reduction potential of the solution.

d. Standard ferrous ammonium sulfate titrant (FAS), approximately 0.10M: Dissolve 39.2 g $\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}$ in distilled water. Add 20 mL conc H_2SO_4 , cool, and dilute to 1000 mL. Standardize solution daily against standard $\text{K}_2\text{Cr}_2\text{O}_7$ digestion solution as follows: Pipet 5.00 mL digestion solution into a small beaker. Add 10 mL reagent water to substitute for sample. Cool to room temperature. Add 1 to 2 drops diluted ferroun indicator and titrate with FAS titrant. Molarity of FAS solution = $[\text{V } \text{K}_2\text{Cr}_2\text{O}_7 \times 0.1] / (\text{VFAS})$ (3) Where: V $\text{K}_2\text{Cr}_2\text{O}_7$ = volume of $\text{K}_2\text{Cr}_2\text{O}_7$ (mL); VFAS = volume of FAS (mL)

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Steps:

1. Wash culture tubes and caps with 20% H_2SO_4 before using to prevent contamination.
2. Place sample in culture tube or ampule and add digestion solution. Carefully run sulfuric acid reagent down inside of vessel so an acid layer is formed under the sample-digestion solution layer and tightly cap tubes or seal ampules, and invert each several times to mix completely.
3. Place tubes or ampules in block digester preheated to 150°C and reflux for 2 h behind a protective shield. CAUTION: These sealed vessels may be under pressure from gases generated during digestion. Wear face and hand protection when handling and dangerous pressures will be generated at 150°C .
4. Cool to room temperature and place vessels in test tube rack. Some mercuric sulfate may precipitate out but this will not affect the analysis.
5. Remove culture tube caps and add small TFE-covered magnetic stirring bar. If ampules are used, transfer contents to a larger container for titrating.
6. Add 0.05 to 0.10 mL (1 to 2 drops) ferroin indicator and stir rapidly on magnetic stirrer while titrating with standardized 0.10M FAS. The end point is a sharp color change from blue-green to reddish brown, although the bluegreen may reappear within minutes. In the same manner reflux and titrate a blank containing the reagents and a volume of distilled water equal to that of the sample.
7. COD is given by

$$\text{COD as mg/L O}_2\text{/L} = [(A-B) \times M \times 8000] / (V_{\text{sample}}) \quad (4)$$

Where: A = volume of FAS used for blank (mL);

B= volume of FAS used for sample (mL);

M=molarity of FAS; 8000= miliquivalent weight of oxygen $\times 1000$ mL/L

Analyze samples in duplicates because of small sample size.

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